

15 - 20 mmHg* Knee High Stockings

Helps energize tired legs and feet, improve circulation and prevent swelling. Soft comfort band is non-restrictive while staying in place all day. Reinforced toe and heel for durability. Sheer, two-way stretch nylon and spandex.



Size	Black	Sand	Natural
Small	68128	68200	68193
Medium	68124	68201	68194
Large	68125	68202	68195
X-Large	68132	68203	68196
XX-Large	68121	68204	68197

Ankle, calf and length are required for proper fit.

Size	Ankle	Calf	Length
Small	6.25"-7.25" (16-18cm)	11"-13.50" (28-34cm)	12"-18" (31-46cm)
Medium	7.50"-9" (19-22cm)	13"-15.50" (33-39cm)	13"-19" (33-48cm)
Large	9.25"-10.75" (24-27cm)	15"-17.50" (38-44cm)	14"-20" (36-51cm)
X-Large	11"-11.50" (28-29cm)	17.75"-18.50" (45-47cm)	15"-22" (38-56cm)
XX-Large	11.75"-12.75" (30-32cm)	18"-19.50" (46-50cm)	15"-22" (38-56cm)



Proper Sizing Directions

Ankle, calf and length are required for proper fit.

- Measure the circumference of the ankle at the narrowest point directly above the ankle bone.
- Measure the circumference of the calf at the widest point.
- Measure the length of the leg from the floor to approx. 1 inch or 2 fingers below the bend of the knee.
- Shoe size provided as a reference only.

Sizing Tips

- Measurements should be taken when you first get up in the morning before swelling can occur. If measurements need to be taken after you have been on your feet for some time, try to elevate your feet for 15 minutes prior to measuring.
- Measurements should be taken while standing. Try to have someone assist you and take the measurements for you as you stand flat footed.
- Take snug measurements against the skin. Try to avoid measuring over slacks or socks.



- Ankle _____
- Calf _____
- Length _____
- Shoe _____



SURGICAL HOUSE
Healthcare Consumables & Equipment

46 KING EDWARD ROAD, OSBORNE PARK 6017

PO BOX 1537, OSBORNE PARK DC WA 6916

surgical@surgicalhouse.com.au

TEL. +61 8 9381 4199

FAX. +61 8 9382 3009

www.surgicalhouse.com.au