

FLOCARE® ENTERAL FEEDING PUMP

A portable enteral feeding pump, designed for controlling the delivery rate of enteral feeds into the gastrointestinal tract

FEATURES

Flocare® Infinity™ III Pump

- Lightweight (430g)
- Compact (154mm x 120mm x 49mm)
- 24-hour battery life
- 30-day feed history
- Multilingual 18 Languages
- Save up to three personalised therapies
- Three therapy modes; Intermittent, Continuous and Volume/Rate
- Flow rate 1-600ml/h in 1ml increments,
- 3-inch Colour display
- Easy to clean, water resistant
- 1-touch manual priming to prepare the set
- · Unique information button with onscreen advice
- Extensive audio and visual safety alarms
- Continuous memory of feeding program and volume delivered
- Lithium ion battery 3.7V DC 3200 mAh,
 6-hour charging time
- · Keypad can be locked
- Accuracy flow rate of +/ 5% with appropriate Flocare Infinity pump giving set
- · Suitable for both adult and paediatric use



Code: **121825** Units per carton: **1**

Directions for Use

- Service check is required at least every 2 years.
- Visit our website to try our Infinity III Pump simulator and view other training resources: www.nutriciaflocare.com/Flocare2020/Infinity3.php

Storage

- Store in a clean, dry environment.
- · Protect from dust, rain and direct sunlight.

Ordering Information

To order contact Nutricia Customer Care 1800 889 480.



Nutricia Australia Pty Ltd The Zenith, Tower B, Level 12, 821 Pacific Hwy, Chatswood NSW 2067 www.nutricia.com.au/adult