

Continence Management Catalogue





MoliCare

HARTMANN has been providing healthcare products for over 200 years and incontinence products for over 50 years. We take a holistic approach to healthcare, placing emphasis on the tangible benefits for patients and those who care for them, aiming to enhance their well-being.


The MoliCare[®] Difference

MoliCare is dedicated to meeting the most stringent standards of quality, dependability, and effectiveness by utilising our advanced technology to develop innovative solutions for managing incontinence. The aim is to address the needs of individuals living with incontinence, enhance the daily duties of skilled caregivers, and satisfy the expectations of healthcare professionals.

Incontinence can disrupt the skin's natural barrier and acid protection mantle, leaving it more susceptible to irritation and damage from urine and faeces.¹ In order to prevent Incontinence Associated Dermatitis (IAD), which can cause skin inflammation, blisters, itchiness or skin erosion²⁻⁴, it is important to take measures to protect the skin from the harmful effects of incontinence.

MoliCare provides absorbent products that are specifically designed to keep the skin dry and maintain its acid protection layer. This can help alleviate the negative consequences of stress factors associated with incontinence and prevent the

development of IAD. This is in line with expert recommendations, which emphasises the importance of preventing IAD as a means of improving quality of life for individuals living with incontinence.⁵



View MoliCare's optimal skin compatibility video to understand how it works.


Scan the QR code or visit: I.ead.me/ph-balance

121212121212



1. Shigeta, Yoshie, et al. "Exploring the relationship between skin property and absorbent pad environment." Journal of Clinical Nursing 18.11 (2009): 1607-1616. 2. Gray, Mikel. "Optimal management of incontinence-associated dermatitis in the elderly." American journal of clinical dermatology 11.3 (2010): 201-210. 3. Gray, Mikel, et al. "Incontinence-associated dermatitis: a consensus." Journal of Wound Ostomy & Continence Nursing 34.1 (2007): 45-54. 4. Junkin, Joan, and Joan Lerner Selekof. "Beyond "diaper rash": Incontinence-associated dermatitis: Does it have you seeing RED?. "Nursing2020 38.11 (2008): 56hn1-56hn10. Warner RR, Stone KJ, Boissy YL. Hydration disrupts human stratum corneum ultrastructure. J Invest Dermatol. 2003;120(2):275–84. 5. Gray, Mikel, et al. "Incontinence-associated dermatitis: a comprehensive review and update." Journal of Wound Ostomy & Continence Nursing 39.1 (2012): 61-74. 6. Bliss et al., An economic evaluation of four skin damage prevention; J Wound Ostomy Continence Nurse. 2007. 7. Data on file. The protective effect of the new diaper design on the epidermal barrier. 2022.

pH of Skin on Incontinence Briefs⁶



MoliCare Absorbent Products:



Significantly lower skin pH values compared to standard briefs, when skin is exposed to alkaline solution.^{6,7}




Lower the prevalence of erythema^{6,7} which is skin redness, an early sign of skin irritation.

Graph data from Bliss Study Absorbent Briefs Containing Curly Fiber Lower (Acidify) Skin pH Reducing Risk for Incontinence Associated Skin Damage in Older Nursing Home Residents.⁶ ANOVA F=8.83 for thigh and F=10.12 for forearm pH differences, p <.001 for both; in posthoc comparisons, skin pH values after wet brief using MoliCare are lower than all other conditions (p=.001 for all).6



- Facilitates rapid fluid distribution within the brief





The MoliCare[®] Skintegrity System

MoliCare offers a comprehensive approach to incontinence management by combining MoliCare absorbent products with MoliCare Skin products.

This combination, known as the MoliCare Skintegrity System, helps to maintain a healthy skin pH balance and helps prevent the occurrence of Incontinence Associated Dermatitis (IAD).

Health Economics of MoliCare

Incidence

- In an acute care study of 5,432 patients¹:
- 45.7% of patients had a prevalence of IAD
- 73% of IAD patients developed IAD following admission

IAD accompanied by unrelieved pressure results in a **5-fold increase** in pressure injury risk²

Cost - IAD¹

The average cost to treat one mild to moderate IAD episode in one resident is \$220.

IAD occurrence can be up to 52%, resulting in costs of \$3,294 per month.

(This example assumes 20% in a 100 bed long term care facility with 75% incontinence at April 2018 rate¹)

Cost - Pressure Injuries^{2,3}





*Except MoliCare® Premium Men Pad 3 Drops, MoliCare® Skin Body Lotion, MoliCare® Hand Cream, MoliCare® Protection Foam

MoliCare absorbent products work in conjunction with MoliCare Skin care products to maintain the acid protection mantle and provide a skin-neutral pH balance.*

Our goal is to provide innovative solutions that improve the management of incontinence, allowing individuals to live with greater comfort and dignity.



3 easy steps to select the right continence aid





STYLE

Multiple styles to meet diverse needs

Select a style based on the level of mobility.

Active

Style based on personal preference.



MoliCare Premium Pads and Pants

Active with restricted mobility and/or cognitive difficulties (e.g. Dementia)

Pull-up recommended looks and feels like regular underwear.

Restricted/supported mobility or bedbound

All-in-one product recommended - easiest option for carers.



MoliCare Premium Mobile



MoliCare Premium Elastic, Form and Fixpants

ABSORBENCY

Stay dry with the right absorption level

- Select an absorbency that matches \bigcirc the level of incontinence.
- Look for the drop rating on the pack. $\langle \mathbf{V} \rangle$

Less Absorbent



\bigotimes	Se pa
\bigotimes	Le too







For optimal comfort and leakage protection

Measure around the widest part between waist and hips.



elect the product size based on the ack's size guide.

eakage may occur if the product is o large or too small.

MoliCare

MoliCare Assortment Overview

		PADS		SUPPORT PAD		PULL-UPS		LARGE PAD		ALL-IN-ONES		BED MATS
	MoliCare Premium Lady Pad	MoliCare Premium Men Pad	MoliCare Pad	MoliCare Rectangular	MoliCare Premium Lady Pants	MoliCare Premium Men Pants	MoliCare Premium Mobile	MoliCare Premium Form	MoliCare Premium Elastic	MoliCare Slip XS	HARTMANN Dignity Bariatric	MoliCare Premium Bed Mat
	Ð	G	Ð		9	0				S.	The second secon	\square
Gender and Mobility Level Absorbency Level	ÔŔ	ŶÂ	ŶŶŔ	够ጰ፟፟፟፟፟ጱ፟ዀ፞ዀ፟҈ዀ	Å 🕅 Å Å h	ំ	0 0 x ki	₿Ŷ\$Ŕ®	ĴĴ Å 🍋	₽9ÅÅ	ĴĴ & 🍋	够ጰ፟፟ዂ፞ዀ፞ዀ
1	Z											
2												
3		D										
4		D						J				
5		D			M C	M L	S M L XL	J				40 x 60cm 60 x 90cm
6							XS M L XL	J	S L XL	xs 👽		
7					M	M C		J	S M L XL			60 x 90cm
8							S M L XL	J ‡	S L XL			
9								J	S L XL			
10							L		S L XL			
Special range								STOOL			2	

	FIXP	ANTS
	MoliCare Premium Fixpants Long Leg	MoliCare Premium Fixpants Short Leg
	Num	
SIZE	COL	OUR
S	YEL	LOW
М	BL	UE
L	BRC	NMO
XL	GRI	EEN
XXL	GR	EY
XXXL	ORANGE	
5XL	PURPLE	

Category symbols

ĵζ

Unisex products

XS - 5XL Available sizes



Products for wearers with restricted mobility Products specifically Products specifically for female wearers

0 Ŵ



absorbency levels



Products for wearers with supported mobility ô

for male wearers



Products specifically for active wearers

Products for bed bound wearers



Due to different cuts, shapes and sizes of the absorbent core in each product group, the exact absorption level differs within each drop level, (i.e. 1 drop does not equal differs within each drop level, (i.e. 1 drop does not equal a defined ISO absorbency). For the exact absorbency measures, refer to specific product information throughout the catalogue. Colours correspond to colour coding of packaging.

MoliCare[®] Fitting Guide Find the right fit



Due to different cuts, shapes and fixation systems, different product categories follow different sizing guides. Measure waist at the largest width between the waist and hips and select the product size based on the size guide printed on the product packaging or below.

Product Group	Absorbency	cm Size	35 4	0 45	50	55 6	60 65	70	75	80 8	85 9	90 95	100	105 110	115	120	125 130	0 135	140	145	150 15	55 160	165	170	175	180 ⁻	185 1	90 19	5 200	200+
MoliCare Premium Lady	5 7	м										M (80	- 120cm	ı)																
& Men Pants		L													L (100 - 1	150cm)													
	6	xs			XS (4	5 - 70c	cm)																							
	568	S					\$	6 (60 -	90сп	n)																				
MoliCare Premium Mobile		м										M (80	- 120cm	J)																
	5680	L													L ('	100 - 1	150cm)													
		XL																	XL	(130 -	170cm)									
		s						S	(70 -	90cm)																				
MoliCare	678	м										M (8	35 - 120	cm)																
Premium Elastic	9 10	L														L (1	15 - 145	5cm)												
		XL																		х	L (140	- 175cn	n)							
MoliCare Slip XS	6	xs		XS (40	- 60сп	n)																								
HARTMANN Dignity Bariatric		-																							Bariat	ric (16	60 - 23	39cm)		

MoliCare Pre	emium Fixpants	cm Size	35	40	45	50	55	60	65	70	75	80	85	90	95	100	105	110	115	120	125	130	135	140	145	150	155	160	165	170	175	180	185	190	195	200	200+
Long Leg	Short Leg																																				
Ye	ellow	s					S (4	15 - 80	Ocm)																												
В	Blue	м								М	(60 -	100c	:m)																								
Br	rown	L												L	(80 -	120ci	m)																				
Gr	reen	XL																		XL	(100	- 160	cm)														
G	Brey	XXL																								XXL	(140	- 180	cm)								
Orange	-	XXXL																												XXX	L (160) - 20	Ocm)				
Purple	-	5XL																																			5XL • 240cm)

PADS

1-5 08

MoliC () () () () () () () () () ()		emium L	ady Pad		Pice Vers
1 Drop	Code	Fit	ISO Max Absorbency in mL*	Units per Pack	Pack per Carton
	168132	Women	181	14	12
2 Drops					
MoliCare'	168634	Women	331	14	18
3 Drops					
MolErre"	168644	Women	504	14	12
4 Drops					
Mol Care"	168682	Women	782	14	12
5 Drops					
MoliCare	168670	Women	1029	14	12

PADS



2 Drops	Code	Fit	Ab
Mel	168600	Men	
3 Drops			
MEN	168603	Men	
4 Drops			
MARCAN	168705	Men	
5 Drops			
	168801	Men	



* Measurement in accordance with ISO 11948/1 (Rothwell method)

MoliCare® Premium Lady Pad 1 Drop

MoliCare® Premium Lady Pad 2, 3, 4, 5 Drops



PADS

2-4 00%

MoliCare Pad 3

2 Drops	Code	Fit	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
Pad D	168100	Unisex	321	28	12
3 Drops					
Matcar" Pad	168101	Unisex	481	28	6
4 Drops					
MatiCarr ^a	168102	Unisex	865	28	6

♣ ÎÅÂÂÂ ⊨ **SUPPORT PADS** MoliCare Rectangular ISO Max Units Pack 3.5 Drops Code Dimensions Absorbency in mL per Pack per Carton 691 56 148530 11 x 35 cm 4


PULL-UPS

MoliCare Premium Lady Pants



Car	915868	L	100 - 150 cm	844	/	8
7 Drops						
ManCover	915877	М	80 - 120 cm	1442	8	4
Ter	915870	L	100 - 150 cm	1373	7	4
A Columb					7 Available	e in Outpatient market

MoliCare Premium Men Pants



១ ៣-៤ ពុំសំភ្លំភ្នំា

nent	ISO Max Absorbency in mL	Units per Pack	Pack per Carton	
cm	843	8	8	
cm	844	7	8	



nent	ISO Max Absorbency in mL	Units per Pack	Pack per Carton	
cm	974	8	4	
cm	1001	7	4	

cm	1455	8	4	
cm	1412	7	4	

Available in Outpatient markets

PULL-UPS

5 – ាំ xs - xL ពុំខ្ញុំ ភ្នំភ្នំរ

MoliCare Premium Mobile





5 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
100	915861	S	60 - 90 cm	941	14	4
Mol Care	915852	М	80 - 120 cm	1017	14	3
	915853	L	100 - 150 cm	1198	14	4
	915854	XL	130 - 170 cm	1183	14	4

6 Drops						
1000	915840	XS	45 - 70 cm	1361	14	4
Moli Care	915831	S	60 - 90 cm	1475	14	4
Martin -	915832	М	80 - 120 cm	1662	14	3
	915833	L	100 - 150 cm	1963	14	4
	915834	XL	130 - 170 cm	2140	14	4

8 Drops						
and the second s	915871	S	60 - 90 cm	1791	14	4
Moster -	915872	М	80 - 120 cm	2015	14	3
	915873	L	100 - 150 cm	2279	14	4
	915874	XL	130 - 170 cm	2421	14	4

10 Drops						
_	915878	М	80 - 120 cm	2626	14	3
P MostCare	915879	L	100 - 150 cm	2613	14	4
	915880	XL	130 - 170 cm	2757	14	4

LARGE PADS

MoliCare Premium Form



4 Drops	Code	Fit	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
MoliCare	168404	Unisex	1493	32	4
5 Drops					
MoliCare'	168405	Unisex	1662	32	4
6 Drops					
MolCare	168406	Unisex	2353	32	4
6 Drops					
O MolCare	168401	Men	2441	32	4
7 Drops					
MoliCare	168407	Unisex	2598	32	4
8 Drops					
MoliCare	168408	Unisex	3017	32	4
9 Drops		Size			
Mailar	168409	Unisex	3497	16	4

4-9 00 x k k k



LARGE PADS

MoliCare®

MoliCare Premium Form + SIZE



	_	
tested	Astra protect	



8 Drops	Code	Fit	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
	168519	Unisex	3347	18	3
Form O AP					

FIXPANTS

MoliCare Premium Fixpants Long Leg



	Code	Size	Waist Measurement	Colour	Units per Pack	Pack per Carton
	947790	S	45 - 80 cm	Yellow	25	1
MoliCare	947791	М	60 - 100 cm	Blue	25	1
(212-212)	947792	L	80 - 120 cm	Brown	25	1
BUG-Gr	947793	XL	100 - 160 cm	Green	25	1
	947794	XXL	140 - 180 cm	Grey	25	1
	947789	XXXL	160 - 200 cm	Orange	25	1
	947825	5XL	200 - 240 cm	Purple	25	1

MoliCare Premium Fixpants Short Leg



	Code	Size	Waist Measurement	Colour	Units per Pack	Pack per Carton
	947711	S	45 - 80 cm	Yellow	25	1
MoliCare	947712	М	60 - 100 cm	Blue	25	1
	947713	L	80 - 120 cm	Brown	25	1
	947714	XL	100 - 160 cm	Green	25	1
	947715	XXL	140 - 180 cm	Grey	25	1



168197 Unisex 1301 32 4	STOOL	Code	Fit	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
	-	168197	Unisex	1301	32	4



With MoliCare®





ALL-IN-ONE

ြ – ြာ S-XL ဂိုဂို မိုမိုေ

MoliCare Premium Elastic





6 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
	165271	S	70 - 90 cm	1667	30	3
MoliCare Month	165272	М	85 - 120 cm	2159	30	3
	165273	L	115 - 145 cm	2349	30	3
	165274	XL	140 - 175 cm	2786	14	4

7 Drops						
	165371	S	70 - 90 cm	1998	30	3
MoliCare Energy	165372	М	85 - 120 cm	2414	30	3
	165373	L	115 - 145 cm	2719	30	3
	165374	XL	140 - 175 cm	3180	14	4

8 Drops						
	165471	S	70 - 90 cm	2356	26	3
MoliCare	165472	М	85 - 120 cm	3144	26	3
	165473	L	115 - 145 cm	3220	24	3
	165474	XL	140 - 175 cm	3591	14	4

9 Drops						
	165571	S	70 - 90 cm	2546	26	3
MoliCare	165572	Μ	85 - 120 cm	3521	26	3
······································	165573	L	115 - 145 cm	3856	24	3
	165574	XL	140 - 175 cm	3951	14	4

10 Drops						
	165671	S	70 - 90 cm	2771	22	3
MoliCare	165672	Μ	85 - 120 cm	3696	14	4
	165673	L	115 - 145 cm	4279	14	4
	165674	XL	140 - 175 cm	4380	14	4

ALL-IN-ONE

MoliCare Slip XS



6 Drops	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units per pack	Pack per carton
	169190	XS	40 - 70 cm	1367	30	4
MoliCare						

HARTMANN Dignity Bariatric



Bariatric	Code	Size	Waist Measurement	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
Dignity Bariatric Brief	916680	Bariatric	160 - 239 cm	2445	8	4
And the second second second						
9						
CO						

6 XS & BARIATRIC η̈́Å Å 🚈





BED MATS

🗴 🗘 ប៉ុនំភំភំាំំំំ 🚝

Notes

MoliCare Premium Bed Mats





5 Drops	Code	Dimensions	ISO Max Absorbency in mL	Units per Pack	Pack per Carton
	161061	40 x 60 cm	449	30	6
Moli Care*	161065	60 x 90 cm	971	30	4
Bod Mat					



Bed Mat

7 Drops					
	161070	60 x 90 cm	1241	25	4







Order Samples


Scan the QR code or visit: http://l.ead.me/molicare-samples



ÔŔ

MoliCare Premium Lady Pad						
Absorbency	Code	Size	Units Per Sample Pack			
	168430	-	2			
2 Drops	168433	-	2			
3 Drops	168434	-	2			
4 Drops	168437	-	2			
5 Drops	168438	-	2			

Mattan	MoliCare Premium Lady Pants				
(Lav	Absorbency	Code	Size	Units Per Sample Pack	
	5 Drops	915810	М	1	
all some		915811	L	1	
0 0 0		915815	М	1	
Ô X ĥ1	7 Drops	915816	L	1	
		7 Available	in Outr	natient markets	





ô	Ř	
~		

Absorbency	Code	Size	Units Per Sample Pack
2 Drops	168503	-	2
3 Drops	168604	-	2
4 Drops	168501	-	2
5 Drops	168609	-	2



MoliCare Pad					
Absorbency	Code	Size	Units Per Sample Pack		
2 Drops	168728	-	2		
3 Drops	168828	-	2		
4 Drops	168928	-	2		

MEN	
បិភិវិតំា	

ំសំភំភំរំា

MoliCare Premium Men Pants					
Absorbency	Code	Size	Units Per Sample Pack		
	915802	М	1		
5 Drops	915803	L	1		
	915807	М	1		
7 Drops	915808	L	1		
— • • • • • • • • • • • • • • • •					

Available in Outpatient markets

MoliCare Premium Mobile					
	Absorbency	Code	Size	Units Per Sample Pack	
		915615	S	2	
		915616	М	2	
191	5 Drops	915617	L	2	
		915648	XL	2	
		915618	XS	2	
		915619	S	2	
	6 Drops	915620	М	2	
		915621	L	2	
		915622	XL	2	
		915623	S	2	
	8 Drops	915624	М	2	
		915625	L	2	
		915656	XL	2	
		915825	М	2	
	10 Drops	915826	L	2	
		915829	XL	2	



Moli	Care	P	rem	iur	n	Fo	٢٢	n

Absorbency	Code Size		Units Per Sample Pack
4 Drops	168204	-	2
5 Drops	168205	-	2
6 Drops	168206	-	2
G Drops O for Men	168201	-	2
7 Drops	168207	-	2
8 Drops	168208	-	2
9 Drops	168209	-	2
+ SIZE	168974	-	2
STOOL	168198	-	1



MoliCare Premium Fixpants Long Leg

<u>=====</u>					
Code	Size	Units Per Sample Pack			
Small	947800	1			
Medium	947801	1			
Large	947802	1			
XL	947803	1			
XXL	947804	1			



MoliCare Premium Elastic				
Absorbency	Code	Size	Units Per Sample Pack	
	165054	S	1	
	165055	М	1	
6 Drops	165056	L	1	
	165057	XL	1	
	165058	S	1	
	165090	М	1	
7 Drops	165059	L	1	
	165060	XL	1	
	165061	S	1	
	165062	М	1	
8 Drops	165063	L	1	
	165064	XL	1	
	165065	S	1	
A 5	165066	М	1	
9 Drops	165067	L	1	
	165068	XL	1	
	165069	S	1	
	165070	М	1	
10 Drops	165071	L	1	
	165072	XL	1	



MoliCare Slip XS

Absorbency	Code	Size	Units Per Sample Pack
6 Drops	169160	XS	2

MoliCare[®] Advanced Features Appendix



Anatomical fit

extra for men

it most.

Discreet

underwear feeling

fixation pants made of

high-quality soft and

Underwear-like

close material

 (\mathbf{v})

Feels dry &

comfortable

Up to 64% faster urine

Extra protection in the

front where men need



Air active

Active Skin

skin protection

Protection features

for optimum skin

wearing comfort.

conditions and

Elasticated

side panels

panels in the

Elasticised side

hip area ensure

to body shape and movements,

delivering

better adaptation

improved comfort for the user.

Hypoallergenic &

dermatologically

The materials used

have extremely low

potential to trigger

allergic reactions.

Secure fixation

Designed to be

underwear.

worn with regular

tested

Comfortable fit High wearing comfort and ideal fit to body.



Ergonomic and easy use **Approximately 20%** faster application as well as easier and more ergonomic handling for caregivers.



Hygienic care Single-use secure protection.



Secure fixation Secure fixation of anatomically shaped incontinence pads.



Wetness indicator on the back sheet Shows when it is time to change.



Active skin protection Skin-friendly, pH balanced ensures maximum skin compatibility. Keeps skin dry with a 3-layer absorbent core and

advanced fast dry effect.



Adjusted fit Safety and comfort through optimal fit; elasticised side panels ensure better adaptation to body shape and


movements, delivering

improved comfort.

26



Advanced skin protection Thanks to a balanced pH



٩

system

Anti-leakage

Fast absorption

with anti-leakage

system and cuffs



odour

tralisa

Odour neutralisation

Absorption system Providing extra protection.



All-round protection Provided by fast absorption, anti-leakage system and odour neutraliser.



for extra security.

1. W. Schneider, H. Mutschler, R. Haller, C. le Mathe: Test for ergonomic evaluation of incontinence products, Oct 5-8, 2015 Dijon. The reduction from 4 to 2 tapes (one for each side) combined with no side panels in the front simplifies the application by reducing the steps from 8 to 4 compared to a traditional slip with 4 side panels and 4 tapes. Additionally, due to the elasticated side panels, there is less need of readjusting the tapes.



Dermatologically tested Helps to maintain healthy skin.





Re-closable tapes 2 re-closable tapes (one for each side) for easy and secure closing.







Tumble dry at low heat



Dermatologically tested Helps to maintain healthy skin.



Discreet

Continence underwear with a low waistline, designed to fit and feel like regular underwear.

MEN with lower waistline and grey tint, fits and feels like regular underwear.



Extra protection Absorbs liquid and provides extra protection for various surfaces.



Extra odour neutralisation Superabsorbent polymers neutralise unpleasant odours.



Prefolded & compressed For easy use and storage.



Quick dry system Absorbent core holds liquid quickly for a dry feeling.



Soft Soft textile-like back sheet provides comfortable feeling.



Soft & skin friendly



3-layer absorbent core Improved absorption core technology that helps to prevent skin irritations.



50x washable Washable up to 50 times at 60°C in a laundry net.



Helps. Cares. Protects.

For more information visit our website or contact us on 1800 805 839 or au_enquiries@hartmann.info

www.molicare.au

PAUL HARTMANN Pty Ltd. ABN 35 000 099 589 Level 5, 1 Thomas Holt Drive, Macquarie Park NSW 2113 ® Registered trademark PAUL HARTMANN AG

9166444_082024

