

## Technical Specification

### Benefits

- Available in two sizes
- Tilt in space
- Fully adjustable
- Compact & easy to fold
- Lightweight
- Pelvic belt fitted as standard
- Height adjustable ergonomic pram handle
- Shopping basket fitted as standard
- Full range of accessories available
- Crash tested to ISO7176-19
- Quick release rear wheels
- Solid PU tyres

### Accessories

- Leg abductor
- Lateral thoracic supports
- Occipital headrest
- Backrest extension
- Foot Straps
- Grab bar
- Contour tray
- Harness
- Sun & rain cover

Prescription form can be downloaded from our website, [www.ineedawheelchair.co.uk](http://www.ineedawheelchair.co.uk)

PART	SIZE 1		SIZE 2	
	MIN	MAX	MIN	MAX
SEAT WIDTH (pelvic supports fitted)	200mm	270mm	260mm	320mm
SEAT DEPTH	230mm	280mm	290mm	350mm
SEAT TO FLOOR HEIGHT (upright position)	-	490mm	-	520mm
BACKREST HEIGHT (from seat base)	-	510mm	-	570mm
BACKREST HEIGHT (from seat base) with back extender	-	625mm	-	720mm
THORACIC WIDTH	120mm	210mm	170mm	240mm
THORACIC HEIGHT (from seat base)	150mm	340mm	160mm	430mm
FOOTPLATE HEIGHT	130mm	280mm	125mm	330mm
FOOTPLATE DEPTH	-	160mm	-	180mm
PUSH HANDLE HEIGHT (from floor)	910mm	1250mm	925mm	1270mm
TILT IN SPACE ANGLE	15°	50°	15°	50°
OVERALL WIDTH	-	580mm	-	640mm
OVERALL DEPTH	-	790mm	-	820mm
OVERALL BACKREST HEIGHT WITH HEADREST (from floor)	-	1070mm	-	1160mm
WEIGHT	14Kg*	-	16Kg*	-
FOLDED SIZE	800x590x500mm		830x640x500mm	
*Rear Wheels Removed				

### Transportation

The Marley has passed the ISO7176-19 test and is approved for motor vehicle use with a maximum user weight of 32kg.

### Upholstery

The Marley comes in black/grey upholstery which is easily removable and machine washable at 40°C



 marley



Manufactured by:

**RMMIS**



British design and quality



The eye-catching looks of the Marley combines simplicity and effectiveness in one design.

The Marley has been developed to meet the child's individual postural needs, with a practical design for daily family living. A versatile buggy available in two sizes with a full range of accessories that complement the many standard features.

#### Tilt in Space

The Marley has a continuous tilt in space range from 15° to 50° to meet individual's requirements

#### Shopping Basket

Under seat basket for additional storage

#### Foldable

Quickly and easily folded in one simple operation to a compact size ideal for both transportation and storage. Removeable rear wheels

#### Ergonomic Pram Handle

A height adjustable handle that can be angled and easily folded to a compact position

#### Lateral Thoracic Support (optional extra)

For additional trunk support easily adjustable in height and width

#### Seat base and Pelvic Supports

Adjustable in depth and width

#### Foot Supports

Height adjustable and flip up to allow for easy transfer

