

VITAL SERIES


A range of red and white contrasting tableware designed to stimulate the appetite and help people more easily see their food or drink due to the colour contrasting rims.

Each piece is designed for maximum functionality and ease of use, making it easy for anyone with cognitive impairment or other mobility limitations to dine with dignity and independence.



Vital Plate and Bowl

- sloping inner edge helps food be easily scooped
- inconspicuous overhang for added grip and stability
- anti-slip ring on base to stop unwanted movement
- ideal for single handed use
- plate edge higher than a standard plate
- plates suited to main meals and soups
- bowl suited to entrees and desserts, or for smaller portions




CROSS SECTION

Vital Ergo Mug

- conical shaped inner cup
- two broad handles for easy grip and added stability
- shape allows contents to be emptied without tilting head
- ideal for drinking in bed or laying down
- thermodynamic function helps contents stay warm or cool

IMPEDIMENTS WITH CONVENTIONAL MUGS



VITAL ERGO MUG



PLATE 26 CM



PLATE 20 CM



BOWL 15 CM



ERGO MUG

	PLATE 26 CM	PLATE 20 CM	BOWL 15 CM	ERGO MUG
CODE	10521	10144	9657	12489
CAPACITY	760 ml	360 ml	185 ml	200 ml
WEIGHT	495 g	300 g	155 g	147 g
DIMENSIONS	260 x 40 mm	200 x 15 mm	155 x 35 mm	80.5 x 108.5 mm