


# MiniONE<sup>®</sup>

Non-Balloon and Capsule Non-Balloon Button



**Innovative  
Technology**

**Improve  
Comfort**



Innovating. Educating. Changing Lives.

# MiniONE®

## Capsule Non-Balloon Button

### Innovative Technology

The AMT MiniONE® Capsule Non-Balloon Button combines the advantages of a low-profile device with the longer life of a non-balloon device.

- No balloon to fail, so replacements occur less often
- Capsule material is **self-lubricating** which facilitates ease and comfort on insertion
- Pull-away tether releases capsule to easily deploy internal bolster
- **Glow Green** feeding port and feed set connectors glow in the dark


### Improve Comfort

Maximizing comfort is key when it comes to patient, caregiver and clinician satisfaction with an enteral feeding device.

- Capsule design dramatically reduces pain during insertion
- Button shaft is comparable to that of a balloon button
- Easy-to-use ratcheted removal tool distends bolster to small profile for ease of removal
- Internal bolster is smaller than a balloon device (less likely to block the pylorus)

Order #	French Size	Length (cm)	Qty/Box
A2-M1-21410-I	14	1.0	1
A2-M1-21412-I	14	1.2	1
A2-M1-21415-I	14	1.5	1
A2-M1-21417-I	14	1.7	1
A2-M1-21420-I	14	2.0	1
A2-M1-21425-I	14	2.5	1
A2-M1-21430-I	14	3.0	1
A2-M1-21812-I	18	1.2	1
A2-M1-21815-I	18	1.5	1
A2-M1-21817-I	18	1.7	1
A2-M1-21820-I	18	2.0	1
A2-M1-21824-I	18	2.4	1
A2-M1-21830-I	18	3.0	1
A2-M1-21834-I	18	3.4	1
A2-M1-21844-I	18	4.4	1
A2-M1-22017-I	20	1.7	1
A2-M1-22024-I	20	2.4	1
A2-M1-22034-I	20	3.4	1
A2-M1-22044-I	20	4.4	1
A2-M1-22417-I	24	1.7	1
A2-M1-22424-I	24	2.4	1
A2-M1-22434-I	24	3.4	1
A2-M1-22444-I	24	4.4	1

Kit includes removal tool



MiniONE® Capsule Non-Balloon

# MiniONE®

## Non-Balloon Button



### Innovative Technology

The AMT MiniONE® Non-Balloon Button combines the security of a non-balloon device with the enhanced features of a low-profile design.

- Cosmetically appealing as the device is flush against the skin
- Ideal for patients who do not tolerate balloon devices – there is no balloon to fail
- Internal retention bolster provides more holes/openings for better feeding and venting of the stomach
- **Glow Green** feeding port and feed set connectors glow in the dark

### Improve Comfort

Maximizing comfort is key when it comes to patient, caregiver and clinician satisfaction with an enteral feeding device. Discomfort is often the product of some degree of trauma.

- Soft, high grade medical silicone is more friendly surface and easier to clean around
- No hard or bulky fill valve to irritate or erode skin surface
- Easy-to-use ratcheted removal tool distends bolster to small profile for ease of removal
- Internal bolster is smaller than a balloon device (less likely to block the pylorus)


Order #	French Size	Length (cm)	Qty/Box
A2-M1-11410-I	14	1.0	1
A2-M1-11412-I	14	1.2	1
A2-M1-11415-I	14	1.5	1
A2-M1-11417-I	14	1.7	1
A2-M1-11420-I	14	2.0	1
A2-M1-11425-I	14	2.5	1
A2-M1-11430-I	14	3.0	1
A2-M1-11812-I	18	1.2	1
A2-M1-11815-I	18	1.5	1
A2-M1-11817-I	18	1.7	1
A2-M1-11820-I	18	2.0	1
A2-M1-11824-I	18	2.4	1
A2-M1-11830-I	18	3.0	1
A2-M1-11834-I	18	3.4	1
A2-M1-11844-I	18	4.4	1
A2-M1-12017-I	20	1.7	1
A2-M1-12024-I	20	2.4	1
A2-M1-12034-I	20	3.4	1
A2-M1-12044-I	20	4.4	1
A2-M1-12417-I	24	1.7	1
A2-M1-12424-I	24	2.4	1
A2-M1-12434-I	24	3.4	1
A2-M1-12444-I	24	4.4	1

Kit includes removal tool



#### MiniONE® Non-Balloon

- Ideal for patients who don't tolerate balloons
- Easy to insert with insertion tool
- Long life (no balloon to break)



**F = Tube Diameter**  
**cm = Stoma Length**



Example: **14F x 2.5cm**

(F) French Size  
(cm) Centimeter



# MiniONE®

## Compatible Feed Sets

		Connector Options					
		Right Angle			Straight		
		2"	12"	24"	2"	12"	24"
Adapter Options	 Y-Port	A2-80222-ISOSAF	A2-81255-ISOSAF**	A2-82455-ISOSAF	A2-80212-ISOSAF	A2-81212-ISOSAF	A2-82455-ISOSAF
	 Bolus	A2-80211-ISOSAF	A2-81221-ISOSAF	A2-82421-ISOSAF	A2-80211-ISOSAF	A2-81211-ISOSAF**	A2-82411-ISOSAF

All MiniONE® feed sets are compatible with Mic-Key® devices. Note: Feeding sets are supplied assembled. Choose the best combination of length, connector and adapter type for your use.

\*\* All MiniONE kits above are packaged with one A2-81255-ISOSAF Y-Port feed set and one A2-81211-ISOSAF Bolus feed set.