

Bathmaster deltis™



Bathmaster deltis™



Compact, lightweight and reclining, with a wealth of features to ensure long term peace of mind



One of the lightest reclining bath lifts available

The Bathmaster Deltis is easy to assemble with no technical knowledge required. With built in convenient holding points and the heaviest part weighing just 6.8kg (15lb) transport and storage become effortless. The all plastic construction is modern, easy to clean and free from corrosion. Designed to be strong, stable & durable, it gives confidence to the user and ensures hassle free, long term usage.

Compact frame fits virtually all bath tubs

Occupying the optimum amount of space, the compact frame of the Bathmaster Deltis is suitable for nearly all bath tubs, yet maintains a high level of stability when in use. Large side flaps allow easy transfer in and out of the bath, and optional swivelling seats are available for those requiring additional assistance.

Latest battery technology ensures high performance

The small, lightweight hand controller is simple to use and floats when not in use, enabling easy recovery if dropped into the water. The lithium ion batteries have no memory effect and very low self-discharging, giving improved performance and battery life. The latest software allows gentle charging, further preserving the batteries.

Difficulties getting into and out of the bath are quickly and easily overcome when using the Bathmaster Deltis. Simple to take in and out of the bath, it separates into two parts for easier lifting and carrying. High quality materials and construction offer the utmost in safety and comfort with no servicing required.



Features

- 3 Year Guarantee
- Lightweight Construction
- Compact Frame
- Latest Battery Technology
- Service Free Motor



Bespoke, service free motor gives long term durability

Fully waterproof with no built in circuit board, the motor has been specifically designed for long term protection against failure or damage. The hand controller plugs directly into the top of the motor ensuring there are no leads from the motor that could be damaged during transportation of the unit.



Three year guarantee giving peace of mind

With a three year guarantee against defects in materials or workmanship, the user can be reassured that the bath lift will be in full working order when they need to use it.



Supportive, reclining backrest offers ultimate relaxation

The gently contoured seat and backrest give both comfort and support when lowering and reclining. The user selects their preferred angle of recline between 10 and 40° allowing optimum positioning for their requirements.



Please order from

Bathmaster Deltis Ordering Details:

Each Bathmaster Delta is supplied complete, including bath lift, hand controller, recharger, side flaps and side flap protectors.

- 09 131 6736 Bathmaster Deltis UK Complete
- 09 132 4078 Bathmaster Deltis Europe Complete
- 09 132 4086 Bathmaster Deltis Australasia Complete
- 09 132 4144 Bathmaster Deltis UK Complete With Blue Covers
- 09 132 4151 Bathmaster Deltis Europe Complete With Blue Covers
- 09 132 4169 Bathmaster Deltis Australasia Complete With Blue Covers
- 09 132 4185 Bathmaster Deltis UK Complete With White Covers
- 09 132 4193 Bathmaster Deltis Europe Complete With White Covers
- 09 132 4201 Bathmaster Deltis Australasia Complete With White Covers

Technical Specifications

Seat size	500 x 375mm (19¾ x 14¾") excluding side flaps
	500 x 695mm (19¾ x 27¼") including side flaps
Backrest size	680 x 350 (26¾ x 13¾")
Seat height	69 to 455mm (2¾ to 18")
Footprint	560 x 295mm (22 x 11½")
Seat weight	6.8kg (15lbs)
Backrest weight	3.7kg (8¼lbs)
Maximum user weight	140kg (22st)
Average number of lifts	9 (based on 80kg load)
Guarantee	3 years on all parts

Features

- 3 Year Guarantee
- Lightweight Construction
- Compact Frame
- Latest Battery Technology
- Service Free Motor



Accessories and Spares

- 09 132 4557 Optional Blue Covers
- 09 132 4565 Optional White Covers
- 09 132 4532 Swivel seat
- 09 132 4540 Swivel transfer seat
- 09 132 5018 Hand Controller
- 09 132 5034 UK Recharger
- 09 132 5026 Europe Recharger
- 09 132 5042 Australasia Recharger

