



Bivona™ Silicone Tracheostomy Tubes

Product Catalogue

Contents

Introduction.....	4
Bivona TTS™.....	5
Bivona TTS Adult Tracheostomy Tubes (Cuffed)	5
Bivona TTS Neonatal Tracheostomy Tubes (Cuffed "V" Neck Flange)	5
Bivona TTS Neonatal Tracheostomy Tubes (Cuffed Straight Neck Flange)*	5
Bivona TTS Paediatric Tracheostomy Tubes (Cuffed "V" Neck Flange)*	5
Bivona TTS Paediatric Tracheostomy Tubes (Cuffed Straight Neck Flange)*	6
Bivona TTS FlexTend™.....	7
Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuff V Neck Flange)	7
Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuff Straight Neck Flange)	7
Bivona TTS Paediatric FlexTend Tracheostomy Tubes (Cuff V Neck Flange)	7
Bivona TTS Paediatric FlexTend Tracheostomy Tubes (Cuff Straight Neck Flange)	7
Bivona TTS FlexTend Plus™.....	8
Bivona TTS Paediatric FlexTend Plus Tracheostomy Tubes (Cuff V Neck Flange)	8
Bivona TTS Paediatric FlexTend Plus Tracheostomy Tubes (Cuff Straight Neck Flange)	8
Bivona FlexTend Plus.....	9
Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)	9
Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes (Extra Length, Straight Neck Flange)*	9
Bivona Cuffless Paediatric FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)*	9
Bivona Cuffless Paediatric FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)*	9
Bivona Cuffless FlexTend Paediatric Tracheostomy Tubes (Extra Length, Straight Neck Flange)*	10
Bivona Cuffless FlexTend Paediatric Tracheostomy Tubes (Extra Length, Straight Neck Flange)*	10
Bivona Aire-Cuf™.....	11
Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes.....	11
Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes with Talk Attachment	11
Bivona Aire-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)*	11
Bivona Aire-Cuf Paediatric Tracheostomy Tubes ("V" Neck Flange)*	12
Bivona Aire-Cuf Paediatric Tracheostomy Tubes (Straight Neck Flange)*	12
Bivona Fome-Cuf™.....	13
Bivona Fome-Cuf Adult Tracheostomy Tubes	13
Bivona Fome-Cuf Adult Tracheostomy Tubes (with Talk Attachment)	13
Bivona Fome-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)*	13
Bivona Fome-Cuf Paediatric Tracheostomy Tubes ("V" Neck Flange)*	13

Bivona Adjustable Hyperflex™	14
Bivona TTS Adjustable Neck Flange Hyperflex Tracheostomy Tubes.....	14
Bivona Mid-Range Aire-Cuf Adjustable Neck Flange Hyperflex Tracheostomy Tubes.....	14
Bivona Uncuffed Tube Adjustable Neck Flange Hyperflex Paediatric "V" Neck Flange Tracheostomy Tubes*	14
Bivona Fixed Flange Hyperflex	15
Bivona Cuffless Extra Length Fixed Hyperflex Tracheostomy Tubes	15
Bivona TTS Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes.....	15
Bivona Cuffed Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes.....	15
Bivona Cuffless/Uncuffed	16
Bivona Cuffless Adult Tracheostomy Tubes	16
Bivona Uncuffed Neonatal Tracheostomy Tubes ("V" Neck Flange)	16
Bivona Uncuffed Neonatal Tracheostomy Tubes (Straight Neck Flange)*	16
Bivona Uncuffed Paediatric Tracheostomy Tubes ("V" Neck Flange)*	16
Bivona Uncuffed Paediatric Tracheostomy Tubes (Straight Neck Flange)*	17
Bivona Customised Silicone Tracheostomy Tubes	18
Accessories	18

Each tube comes individually packaged with obturator, disconnection wedge, and twill tape.

Bivona Silicone Tracheostomy Tubes

Silicone offers important advantages over conventional tube materials:

Softness and flexibility

For comfort, ease of use, and greater range of mobility.

Biocompatibility

Hydrophobic and hypoallergenic, tissue-friendly and easy to clean.

Durability

Robust for reliable reprocessing, withstands repeated sanitisation up to ten times (adult) and five times (for neonatal and paediatrics).



Bivona silicone tracheostomy tubes are available in a range of tube configuration options to suit each, age, size, and application requirements. Customised Bivona silicone tubes are available.

For comfort

A variety of tubes designed to aid in patient comfort:

- › allows 15 mm connector to be distal from chin.
- › reduces the risk of device-related pressure ulcers.
- › option to utilise self-expanding cuff to mitigate risks of over inflation.

For design

Integrating unique features to aid patient comfort, ease-of-use and avoidance of contamination.

For useability and safety

Hydrophobic material for easy cleaning, free from DEHP, latex and plasticisers.

For cost-efficiency

Bivona tubes can be cleaned and sanitised for up to ten re-insertions for adult sizes or five times for neonatal and paediatric sizes.

For choice

Combining the softness and flexibility of silicone, with an extensive range of size and shape options for all ages.



Bivona TTS

The TTS cuff design offers the benefits of a cuff with the profile of a cuffless tube. When the silicone cuff is totally deflated it adds no significant dimension to the outer diameter of the tube's shaft, yet it can be partially or completely inflated for varied airway management needs.



Hospital



Home



MR
Conditional



670180

Bivona TTS Adult Tracheostomy Tubes (Cuffed)							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
670150	5.0	5.0	7.4	60.0	100°	9.5	1
670160	6.0	6.0	8.8	70.0	100°	11.0	1
670170	7.0	7.0	10.0	80.0	100°	12.0	1
670175	7.5	7.5	10.5	80.0	100°	12.5	1
670180	8.0	8.0	11.0	88.0	100°	13.0	1
670185	8.5	8.5	11.8	88.0	100°	14.0	1
670190	9.0	9.0	12.3	98.0	100°	14.5	1



67N030

Bivona TTS Neonatal Tracheostomy Tubes (Cuffed "V" Neck Flange)							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67N025	2.5	2.5	4.0	30.0	120°	4.0	1
67N030	3.0	3.0	4.7	32.0	120°	4.7	1
67N035	3.5	3.5	5.3	34.0	120°	5.3	1
67N040	4.0	4.0	6.0	36.0	120°	6.0	1



67SN030

Bivona TTS Neonatal Tracheostomy Tubes (Cuffed Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67SN025	2.5	2.5	4.0	30.0	120°	4.0	1
67SN030	3.0	3.0	4.7	32.0	120°	4.7	1
67SN035	3.5	3.5	5.3	34.0	120°	5.3	1
67SN040	4.0	4.0	6.0	36.0	120°	6.0	1



67P045

Bivona TTS Paediatric Tracheostomy Tubes (Cuffed "V" Neck Flange)							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67P025	2.5	2.5	4.0	38.0	120°	4.0	1
67P030	3.0	3.0	4.7	39.0	120°	4.7	1
67P035	3.5	3.5	5.3	40.0	120°	5.3	1
67P040	4.0	4.0	6.0	41.0	120°	6.0	1
67P045	4.5	4.5	6.7	42.0	120°	6.7	1
67P050	5.0	5.0	7.3	44.0	120°	7.3	1
67P055	5.5	5.5	8.0	46.0	120°	8.0	1

Bivona TTS




67SP045

Bivona TTS Paediatric Tracheostomy Tubes (Cuffed Straight Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67SP025	2.5	2.5	4.0	38.0	120°	4.0	1
67SP030	3.0	3.0	4.7	39.0	120°	4.7	1
67SP035	3.5	3.5	5.3	40.0	120°	5.3	1
67SP040	4.0	4.0	6.0	41.0	120°	6.0	1
67SP045	4.5	4.5	6.7	42.0	120°	6.7	1
67SP050	5.0	5.0	7.3	44.0	120°	7.3	1
67SP055	5.5	5.5	8.0	46.0	120°	8.0	1

Bivona TTS FlexTend

Radiopaque silicone tracheostomy tubes with a wire reinforced extension, proximal to the neck flange. The silicone cuff lays flat against the shaft with deflated. Proximal extension allows 15 mm connector to be distal from the chin. The tube has a clear inflation balloon with a one way valve and integral 15 mm connection. It is manufactured with Bivona SuperSlick™ formulation.



67NFP25

Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuff V Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67NFP25	2.5	2.5	4.0	20.0	30.0	120°	4.0	1
67NFP30	3.0	3.0	4.7	20.0	32.0	120°	4.7	1
67NFP35	3.5	3.5	5.3	20.0	34.0	120°	5.3	1
67NFP40	4.0	4.0	6.0	20.0	36.0	120°	6.0	1



67NFPS25

Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuff Straight Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67NFPS25	2.5	2.5	4.0	20.0	30.0	120°	4.0	1
67NFPS30	3.0	3.0	4.7	20.0	32.0	120°	4.7	1
67NFPS35	3.5	3.5	5.3	20.0	34.0	120°	5.3	1
67NFPS40	4.0	4.0	6.0	20.0	36.0	120°	6.0	1



67PFS25

Bivona TTS Paediatric FlexTend Tracheostomy Tubes (Cuff V Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFS25	2.5	2.5	4.0	20.0	38.0	120°	4.0	1
67PFS30	3.0	3.0	4.7	20.0	39.0	120°	4.7	1
67PFS35	3.5	3.5	5.3	20.0	40.0	120°	5.3	1
67PFS40	4.0	4.0	6.0	20.0	41.0	120°	6.0	1
67PFS45	4.5	4.5	6.7	30.0	42.0	120°	6.7	1
67PFS50	5.0	5.0	7.3	30.0	44.0	120°	7.3	1
67PFS55	5.5	5.5	8.0	30.0	46.0	120°	8.0	1




67PFSS25

Bivona TTS Paediatric FlexTend Tracheostomy Tubes (Cuff Straight Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFSS25	2.5	2.5	4.0	20.0	38.0	120°	4.0	1
67PFSS30	3.0	3.0	4.7	20.0	39.0	120°	4.7	1
67PFSS35	3.5	3.5	5.3	20.0	40.0	120°	5.3	1
67PFSS40	4.0	4.0	6.0	20.0	41.0	120°	6.0	1
67PFSS45	4.5	4.5	6.7	30.0	42.0	120°	6.7	1
67PFSS50	5.0	5.0	7.3	30.0	44.0	120°	7.3	1
67PFSS55	5.5	5.5	8.0	30.0	46.0	120°	8.0	1

Bivona TTS FlexTend Plus

Radiopaque silicone tracheostomy tubes with a wire reinforced extension, proximal to the neck flange. The silicone cuff lays flat against the shaft when deflated. Proximal extension allows 15 mm connector to be distal from the chin. The tube has a clear inflation balloon with a one way valve and integral 15 mm connection. FlexTend Plus tubes also provide additional distal length for patients requiring additional overall length. It is manufactured with Bivona SuperSlick™.



67PFP40

Bivona TTS Paediatric FlexTend Plus Tracheostomy Tubes (Extra Length, Cuff V Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFP40	4.0	4.0	6.0	30.0	44.0	120°	4.0	1
67PFP45	4.5	4.5	6.7	30.0	48.0	120°	4.7	1
67PFP50	5.0	5.0	7.3	30.0	50.0	120°	5.3	1
67PFP55	5.5	5.5	8.0	30.0	52.0	120°	6.0	1
67PFP60	6.0	6.0	8.7	30.0	56.0	120°	6.7	1



Bivona TTS Paediatric FlexTend Plus Tracheostomy Tubes (Extra Length, Cuff Straight Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFPS40	4.0	4.0	6.0	30.0	44.0	120°	4.0	1
67PFPS45	4.5	4.5	6.7	30.0	48.0	120°	4.7	1
67PFPS50	5.0	5.0	7.3	30.0	50.0	120°	5.3	1
67PFPS55	5.5	5.5	8.0	30.0	52.0	120°	6.0	1
67PFPS60	6.0	6.0	8.7	30.0	56.0	120°	6.7	1

Bivona FlexTend

Reinforced proximal and distal sections provide a high level of kink and crush resistance. Extended proximal tube allows 15 mm connector to be distal from the chin.



Hospital



Home



MR
Conditional



60NFP30

Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes ("V" Neck Flange)								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60NFP25	2.5	2.5	4.0	20.0	30.0	120°	N/A	1
60NFP30	3.0	3.0	4.7	20.0	32.0	120°	N/A	1
60NFP35	3.5	3.5	5.3	20.0	34.0	120°	N/A	1
60NFP40	4.0	4.0	6.0	20.0	36.0	120°	N/A	1



60NFP35

Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes (Straight Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60NFP35	3.5	3.5	5.3	20.0	34.0	120°	N/A	1



60PFS40

Bivona Cuffless Paediatric FlexTend Tracheostomy Tubes ("V" Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFS25	2.5	2.5	4.0	20.0	38.0	120°	N/A	1
60PFS30	3.0	3.0	4.7	20.0	39.0	120°	N/A	1
60PFS35	3.5	3.5	5.3	20.0	40.0	120°	N/A	1
60PFS40	4.0	4.0	6.0	20.0	41.0	120°	N/A	1
60PFS45	4.5	4.5	6.7	30.0	42.0	120°	N/A	1
60PFS50	5.0	5.0	7.3	30.0	44.0	120°	N/A	1
60PFS55	5.5	5.5	8.0	30.0	46.0	120°	N/A	1



60PFP45

Bivona Cuffless Pediatric FlexTend Plus Tracheostomy Tubes (Extra Length, "V" Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFP40	4.0	4.0	6.0	30.0	44.0	120°	N/A	1
60PFP45	4.5	4.5	6.7	30.0	48.0	120°	N/A	1
60PFP50	5.0	5.0	7.3	30.0	50.0	120°	N/A	1
60PFP55	5.5	5.5	8.0	30.0	52.0	120°	N/A	1
60PFP60	6.0	6.0	8.7	30.0	56.0	120°	N/A	1

Bivona FlexTend



60PFPS50

Bivona Cuffless FlexTend Plus Paediatric Tracheostomy Tubes (Extra Length, Straight Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFPS40	4.0	4.0	6.0	30.0	44.0	120°	N/A	1
60PFPS50	5.0	5.0	7.3	30.0	50.0	120°	N/A	1
60PFPS60	6.0	6.0	8.7	30.0	56.0	120°	N/A	1
60PFPS55	5.5	5.5	8.0	30.0	52.0	120°	N/A	1
60PFPS60	6.0	6.0	8.7	30.0	56.0	120°	N/A	1



60PFSS25

Bivona Cuffless FlexTend Paediatric Tracheostomy Tubes (Straight Neck Flange)*								
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFSS25	2.5	2.5	4.0	20.0	38.0	120°	N/A	1
60PFSS30	3.0	3.0	4.7	20.0	39.0	120°	N/A	1
60PFSS35	3.5	3.5	5.3	20.0	40.0	120°	N/A	1
60PFSS40	4.0	4.0	6.0	20.0	41.0	120°	N/A	1
60PFSS45	4.5	4.5	6.7	30.0	42.0	120°	N/A	1
60PFSS50	5.0	5.0	7.3	30.0	44.0	120°	N/A	1
60PFSS55	5.5	5.5	8.0	30.0	46.0	120°	N/A	1

Bivona Aire-Cuf

A flexible, adjustable tube with a soft secure cuff. Neonatal and paediatric tubes are reinforced to give a high level of kink and crush resistance. Select adult configurations feature a talk attachment allowing vocalisation for ventilator dependent patients with viable vocal cords.



Hospital



Home



MR Conditional



750160

Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
750150	5.0	5.0	7.4	60.0	100°	17.0	1
750160	6.0	6.0	8.8	70.0	100°	18.0	1
750170	7.0	7.0	10.0	80.0	100°	21.0	1
750180	8.0	8.0	11.0	88.0	100°	25.0	1
750190	9.0	9.0	12.3	98.0	100°	28.0	1
750195	9.5	9.5	13.3	98.0	100°	28.0	1



755180

Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes with Talk Attachment							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
755160-SM	6.0	6.0	8.7	70.0	100°	18.0	1
755170-SM	7.0	7.0	10.0	80.0	100°	21.0	1
755180-SM	8.0	8.0	11.0	88.0	100°	25.0	1
755190-SM	9.0	9.0	12.3	98.0	100°	28.0	1
755195-SM	9.5	9.5	13.3	98.0	100°	28.0	1



65N035

Bivona Aire-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65N025	2.5	2.5	4.0	30.0	120°	12.0	1
65N030	3.0	3.0	4.7	32.0	120°	14.0	1
65N035	3.5	3.5	5.3	34.0	120°	14.0	1
65N040	4.0	4.0	6.0	36.0	120°	18.0	1

Bivona Aire-Cuf



65P030

Bivona Aire-Cuf Paediatric Tracheostomy Tubes ("V" Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65P025	2.5	2.5	4.0	38.0	120°	12.0	1
65P030	3.0	3.0	4.7	39.0	120°	14.0	1
65P035	3.5	3.5	5.3	40.0	120°	14.0	1
65P040	4.0	4.0	6.0	41.0	120°	18.0	1
65P045	4.5	4.5	6.7	42.0	120°	18.0	1
65P050	5.0	5.0	7.3	44.0	120°	20.0	1
65P055	5.5	5.5	8.0	46.0	120°	20.0	1




65SP030

Bivona Aire-Cuf Paediatric Tracheostomy Tubes (Straight Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65SP025	2.5	2.5	4.0	38.0	120°	12.0	1
65SP030	3.0	3.0	4.7	39.0	120°	14.0	1
65SP035	3.5	3.5	5.3	40.0	120°	14.0	1
65SP040	4.0	4.0	6.0	41.0	120°	18.0	1
65SP045	4.5	4.5	6.7	42.0	120°	18.0	1
65SP050	5.0	5.0	7.3	44.0	120°	20.0	1
65SP055	5.5	5.5	8.0	46.0	120°	20.0	1

Bivona Fome-Cuf

The self expanding cuff helps reduce the risk of over or under-inflation associated with traditional cuffs. Select adult configurations feature a talk attachment allowing vocalisation for ventilator dependent patients with viable vocal cords.



850180

Bivona Fome-Cuf Adult Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
850150	5.0	5.0	7.3	60.0	100°	20.0	1
850160	6.0	6.0	8.7	70.0	100°	30.0	1
850170	7.0	7.0	10.0	80.0	100°	32.0	1
850180	8.0	8.0	11.0	88.0	100°	33.0	1
850190	9.0	9.0	12.3	98.0	100°	33.0	1
850195	9.5	9.5	13.3	98.0	100°	33.0	1



855180

Bivona Fome-Cuf Adult Tracheostomy Tubes (with Talk Attachment)							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
855150	5.0	5.0	7.3	60.0	100°	20.0	1
855160	6.0	6.0	8.7	70.0	100°	30.0	1
855170	7.0	7.0	10.0	80.0	100°	32.0	1
855180	8.0	8.0	11.0	88.0	100°	33.0	1
855190	9.0	9.0	12.3	98.0	100°	33.0	1



85N030

Bivona Fome-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
85N030	3.0	3.0	4.7	32.0	N/A	14.0	1
85N035	3.5	3.5	5.3	34.0	N/A	14.0	1
85N040	4.0	4.0	6.0	36.0	N/A	18.0	1



85P035

Bivona Fome-Cuf Paediatric Tracheostomy Tubes ("V" Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
85P025	2.5	2.5	4.0	38.0	N/A	12.0	1
85P030	3.0	3.0	4.7	39.0	N/A	14.0	1
85P035	3.5	3.5	5.3	40.0	N/A	14.0	1
85P040	4.0	4.0	6.0	41.0	N/A	18.0	1
85P045	4.5	4.5	6.7	42.0	N/A	18.0	1
85P050	5.0	5.0	7.3	44.0	N/A	20.0	1
85P055	5.5	5.5	8.0	46.0	N/A	20.0	1

Bivona Adjustable Hyperflex


Adjustable flange tubes allow for immediate modification of flange position to meet a variety of patient needs. It can be adjusted for horizontal and vertical shaft drop to accommodate unusual anatomy or pathology. Flexible wire-reinforced silicone shafts provide improved kink and crush resistance.



Hospital



MR
Conditional



67HA60

Bivona TTS Adjustable Neck Flange Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67HA60	6.0	6.0	9.2	110.0	180°	11.0	1
67HA70	7.0	7.0	10.6	120.0	180°	12.0	1
67HA80	8.0	8.0	11.7	130.0	180°	13.0	1
67HA90	9.0	9.0	12.9	140.0	180°	14.5	1



75HA70

Bivona Mid-Range Aire-Cuf Adjustable Neck Flange Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
75HA60	6.0	6.0	9.2	110.0	180°	18.0	1
75HA70	7.0	7.0	10.6	120.0	180°	21.0	1
75HA80	8.0	8.0	11.7	130.0	180°	25.0	1
75HA90	9.0	9.0	12.9	140.0	180°	28.0	1



60HA35

Bivona Uncuffed Tube Adjustable Neck Flange Hyperflex Paediatric "V" Neck Flange Tracheostomy Tubes*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Maximum Usable Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60HA25	2.5	2.5	4.0	55.0	180°	N/A	1
60HA30	3.0	3.0	4.7	60.0	180°	N/A	1
60HA35	3.5	3.5	5.3	65.0	180°	N/A	1
60HA40	4.0	4.0	6.0	70.0	180°	N/A	1
60HA45	4.5	4.5	6.7	75.0	180°	N/A	1
60HA50	5.0	5.0	7.3	80.0	180°	N/A	1
60HA55	5.5	5.5	8.0	85.0	180°	N/A	1

Bivona Fixed Flange Hyperflex

The fixed flange tube provides a more permanent solution for patients who require extra vertical and horizontal length. Flexible wire-reinforced silicone shafts provide improved kink and crush resistance.



Hospital



Home



MR
Conditional



60AFHXL70

Bivona Cuffless Extra Length Fixed Hyperflex Tracheostomy Tubes*

Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60AFHXL60	6.0	6.0	9.2	100.0	180°	N/A	1
60AFHXL70	7.0	7.0	10.6	110.0	180°	N/A	1
60AFHXL80	8.0	8.0	11.7	120.0	180°	N/A	1
60AFHXL90	9.0	9.0	12.9	130.0	180°	N/A	1



67FHXL70

Bivona TTS Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes*

Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67FHXL60	6.0	6.0	9.2	100.0	180°	11.0	1
67FHXL70	7.0	7.0	10.6	110.0	180°	12.0	1
67FHXL80	8.0	8.0	11.7	120.0	180°	13.0	1
67FHXL90	9.0	9.0	12.9	130.0	180°	14.5	1



75FHXL60

Bivona Cuffed Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes*

Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
75FHXL60	6.0	6.0	9.2	100.0	180°	18.0	1
75FHXL70	7.0	7.0	10.6	110.0	180°	21.0	1
75FHXL80	8.0	8.0	11.7	120.0	180°	25.0	1
75FHXL90	9.0	9.0	12.9	130.0	180°	28.0	1

Bivona Cuffless

When an uncuffed tube is indicated, the all silicone construction of this soft and flexible tube ensures comfort for the patient.



Hospital



Home



MR
Conditional



60A160

Bivona Cuffless Adult Tracheostomy Tubes*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60A150	5.0	5.0	7.4	60.0	N/A	N/A	1
60A160	6.0	6.0	8.8	70.0	N/A	N/A	1
60A170	7.0	7.0	10.0	80.0	N/A	N/A	1
60A180	8.0	8.0	11.0	88.0	N/A	N/A	1
60A190	9.0	9.0	12.3	98.0	N/A	N/A	1

* Products listed above are MR Safe.



60N030

Bivona Uncuffed Neonatal Tracheostomy Tubes ("V" Neck Flange)**							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60N025	2.5	2.5	4.0	30.0	120°	N/A	1
60N030	3.0	3.0	4.7	32.0	120°	N/A	1
60N035	3.5	3.5	5.3	34.0	120°	N/A	1
60N040	4.0	4.0	6.0	36.0	120°	N/A	1



60SN025

Bivona Uncuffed Neonatal Tracheostomy Tubes (Straight Neck Flange)**							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60SN025	2.5	2.5	4.0	30.0	120°	N/A	1
60SN030	3.0	3.0	4.7	32.0	120°	N/A	1
60SN035	3.5	3.5	5.3	34.0	120°	N/A	1
60SN040	4.0	4.0	6.0	36.0	120°	N/A	1



60P040

Bivona Uncuffed Paediatric Tracheostomy Tubes ("V" Neck Flange)**							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60P025	2.5	2.5	4.0	38.0	120°	N/A	1
60P030	3.0	3.0	4.7	39.0	120°	N/A	1
60P035	3.5	3.5	5.3	40.0	120°	N/A	1
60P040	4.0	4.0	6.0	41.0	120°	N/A	1
60P045	4.5	4.5	6.7	42.0	120°	N/A	1
60P050	5.0	5.0	7.3	44.0	120°	N/A	1
60P055	5.5	5.5	8.0	46.0	120°	N/A	1

Bivona Cuffless



Bivona Uncuffed Paediatric Tracheostomy Tubes (Straight Neck Flange)*							
Product Code	Tube Size (mm)	Tube I.D. (mm)	Tube O.D. (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60SP025	2.5	2.5	4.0	38.0	120°	N/A	1
60SP030	3.0	3.0	4.7	39.0	120°	N/A	1
60SP035	3.5	3.5	5.3	40.0	120°	N/A	1
60SP040	4.0	4.0	6.0	41.0	120°	N/A	1
60SP045	4.5	4.5	6.7	42.0	120°	N/A	1
60SP050	5.0	5.0	7.3	44.0	120°	N/A	1
60SP055	5.5	5.5	8.0	46.0	120°	N/A	1

Bivona Customised Silicone Tracheostomy Tubes


Treating a patient with difficult or unusual anatomy can be a challenge; getting that patient a proper fitting tracheostomy tube shouldn't be. Bivona customisation lab technicians can hand-craft a tube of soft, durable silicone to conform to the patient's neck, stoma, and trachea.

Designing a Bivona customised tracheostomy tube can be completed using our on-line templates.

Bivona customised tracheostomy tube delivery times:

- › Standard service (sterile): 3-4 weeks
- › Standard service (non-sterile): 7 business days
- › Express service (non-sterile): 3 business days
- › Next-day service (non-sterile): shipped same day if order is received by 10 a.m. CST
- › Next-day service (non-sterile): shipped same day if order is received by 12 p.m. CST (additional fees apply)

Visit bivonacustom.com for more information or contact us at:
Email: bivonacustomorder@icumed.com



Accessories

Accessories		
Product Code	Description	Units Per Case
CON15S	Side Port Connectors NS	50
REDCAP	15 mm Red Cap Plugs	26
DWG300	Disconnect Wedge 15 mm	25
570016	Thermovent™ T (Sterile)	50
100/580/015	Thermovent 600 HME	20
100/570/022	Thermovent T2 (Sterile)	50
100/575/010	Oxygen Delivery for HME Thermovent T Thermovent O ₂	10
100/582/000	Thermovent 1200 HME, Adult, Straight	20

