

Air Lift Chair



The ultimate in comfort and freedom



Features

- Unique modular design
- Unique air system for pressure relief
- Powered adjustment (individual controls)
- Adjustable backrest and leg rest
- Tilt assistance for standing
- Rear wheels for easy transport
- 4 way stretchable incontinent polyurethane cover (fire retardant)
- Ergonomical push handles

Options

- Blue or Beige covers
- Single or dual actuator system

Technical Data

Measurement	Seat Width	Seat Height	Seat Depth	Overall Depth	Overall Height	Overall Width
mm	600	500	500	1000	1100	800
inch	23	20	20	39	43	32
	Reclining		Foot Rest		SWL	
	45° ~105°		5° ~90°		125KG	

Technical Specifications

