

CA-MI

Surgical Suction Unit

**2 YEAR
WARRANTY**

New ASKIR₃₀



▶ *Compact and portable unit*

▶ *Ideal for General Practice
& home use*

TECHNICAL SPECIFICATIONS

| | |
|--|----------------------------|
| Power Feeding other voltages available upon request | 230V-50hz |
| Max Vacuum - adjustable | 0.80 Bar/ -80kPa/ -600mmHg |
| Max free air flow | 40 l/min |
| Power Consumption | 107 VA |
| Fuse | 1x1.6A 250V |
| Intermittent duty cycle at 35 degrees and 110% operating Voltage | 120 ON/ 60 OFF |
| Weight | 3.60 kg |
| Size | 350mm x 210mm x 180mm |
| Warranty | 2 years |
| Shipping carton | 4 |
| Manufacturing | 100% Italy |

CONFIGURATIONS

With AUTOCLAVE JAR *with overflow valve system, antibacterial and hydrophobic filter, complete set of tubes 6 x 10mm, conical connector for probes.

*with 1000ml autoclavable jar: RE310100/02



SURGICAL HOUSE

166 RAILWAY PARADE, LEEDERVILLE WA 6007

PO BOX 84, LEEDERVILLE WA 6902

surgical@surgicalhouse.com.au

TEL. +61 8 9381 4199

FAX. +61 8 9382 3009

www.surgicalhouse.com.au