CodaCare







Model 225L

Spring Adjustable Overbed Table - Laminated Top

- Sturdy steel frame powder coated light beige
- Spring adjustable with locking hand wheel
- Moisture resistant laminated top, 800 x 400mm (18mm thick)
- Impact resistant PVC matching edging
- 50mm twin wheeled castors (inclusive of 2x lockable castors)
- Height adjustment range 740 to 1050 (mm)
- Safe Working Load: 25kg evenly distributed

Table Top Dimensions: W: 800mm - D: 400mm - H: 18mm

Beech Grey Maple Maple Popular Colours:

You can choose from our existing range or have it resized or custom made to your specification.

Model 225T

Spring Adjustable Overbed Table - Thermoformed Top

- Sturdy steel frame powder coated light beige
- Spring adjustable with locking hand wheel
- Gap free Thermoformed top, 800 x 400mm (25mm thick)
- 5mm raised edge to prevent spillages
- 50mm twin wheeled castors (inclusive of 2x lockable castors)
- Height adjustment range 740 to 1050 (mm)
- Safe Working Load: 25kg evenly distributed

Table Top Dimensions: W: 800mm - D: 400mm - H: 25mm

Popular Colours: Maple Beech Grey

You can choose from our existing range or have it resized or custom made to your specification.





Model 225LOW

Spring Adjustable Overbed Table - Low Height Frame

- Sturdy steel frame powder coated light beige
- Spring adjustable with locking hand wheel
- Moisture resistant laminated top, 800 x 400mm (18mm thick)
- Impact resistant PVC matching edging
- 50mm twin wheeled castors (inclusive of 2x lockable castors)
- Height adjustment range 500 to 750 (mm)
- Safe Working Load: 25kg evenly distributed

Table Top Dimensions: W: 800mm - D: 400mm - H: 18mm

Popular Colours: Maple Beech

You can choose from our existing range or have it resized or custom made to your specification.

Model 240

Spring Adjustable Overbed Table – Split Tilt Top

- Sturdy steel frame powder coated light beige
- Spring adjustable with locking hand wheel
- Moisture resistant split tilt top
- Impact resistant PVC matching edging
- 50mm twin wheeled castors (inclusive of 2x lockable castors)
- Height adjustment range 740 to 1050 (mm) Safe Working Load: 25kg evenly distributed

Table Top Dimensions: W: 600mm - D: 400mm (Tilt) / W: 200mm - D: 400mm (Fix)

Popular Colours: Maple Beech

You can choose from our existing range or have it resized or custom made to your specification.





EVO ONE

Premium Overbed Table

- ✓ One hand ergonomic height adjustment
- ✓ Extremely strong powder coated C-base
- ✓ The table top surface is smooth, non-absorbent and free from gaps
- ✓ Heavy duty 70mm castors (inclusive of 2x lockable castors)
- ✓ Height adjustment range 760 to 1020 (mm)
- ✓ Safe Working Load: 20kg evenly distributed
- ✓ 10 year warranty on High Pressure Laminate (HPL)

Table Top Dimensions: W: 800mm - D: 400mm - H: 20mm

Popular Colours: Maple Beech Grey Beige

You can choose from our existing range or have it resized or custom made to your specification.

EVO ONE LOW

Premium Low Height Overbed Table

- One hand ergonomic height adjustment
- ✓ Extremely strong powder coated C-base
- The table top surface is smooth, non-absorbent and free from gaps
- Heavy duty 70mm castors (inclusive of 2x lockable castors)
- Height adjustment range 840 to 880 (mm)
- ✓ Safe Working Load: 20kg evenly distributed
- 10 year warranty on High Pressure Laminate (HPL)

Table Top Dimensions: W: 900mm - D: 400mm - H: 20mm

Popular Colours: Maple Beech Grey Beige

You can choose from our existing range or have it resized or custom made to your specification.





EVO ONE TILT

Premium Overbed Table with Split Tilt Top

- ✓ One hand ergonomic height adjustment
- ✓ Extremely strong powder coated C-base
- ✓ The table top surface is smooth, non-absorbent and free from gaps
- ✓ Heavy duty 70mm castors (inclusive of 2x lockable castors)
- ✓ Height adjustment range 760 to 1020 (mm)
- ✓ Safe Working Load: 20kg evenly distributed
- √ 10 year warranty on High Pressure Laminate (HPL)

Table Top Dimensions: W: 680mm - D: 400mm (Tilt) / W: 280mm - D: 400mm (Fix)

Popular Colours: Maple Beech Grey Beige

You can choose from our existing range or have it resized or custom made to your specification.

PremiumLift

Multipurpose Overbed Table

- ✓ One hand ergonomic height adjustment
- Extremely strong powder coated C-base
- ✓ The table top surface is smooth, non-absorbent and free from gaps.
- 50mm twin wheeled castors (inclusive of 2x lockable castors)
- Height adjustment range 740 to 1050 (mm)
- ✓ Safe Working Load: 20kg evenly distributed
- ✓ Also available with standard non-tilting table top

Table Top Dimensions: W: 600mm - D: 400mm (Tilt) / W: 200mm - D: 400mm (Fix)

Popular Colours: Maple Beech Grey Beige

You can choose from our existing range or have it resized or custom made to your specification.



CodaCare





