



CU-SPI (SPI - IPAD AED)

EASY-TO-USE SEMI-AUTOMATED EXTERNAL DEFIBRILLATOR (AED) THAT IS SMALL, LIGHT & PORTABLE

KEY FEATURES

- Designed for public use and thus can be used by minimally trained individuals
- Ambient noise detector measures level of background noise and adjusts the volume of the voice prompts accordingly
- Integrated pad storage: Pre-connected pads are stored in a clear compartment on the underside of the unit
- Visual indicator shows at a glance battery life, unit status and pad status

- ✓ Easily switch from adult to child mode without changing pads
- ✓ Auto volume adjustment
- ✓ CPR detection
- ✓ Daily, weekly, monthly self checks

5 YEAR WARRANTY EXCLUDES BATTERY

CU-SP1



PRODUCTS

Product Number	Description
CU-SP1	AED iPad SP1 - Semi Automated
CU-SP1-A-BAG-3010	SPI Carrying case
CU-A1007S	Disposable Defib pads (Adult)
CU-A1102S	Disposable Defib pads (Child)
CU-SPI-A08	Wall cabinet
CU-SPI-A12	Wall bracket
CU-SA1103BB	Disposable battery pack
CU-SPI-A10-SD	SD card
CU-SPT	Trainer model
CU-SPT-OA03-BAG	Carrying case (SP1 Trainer)



Wall Cabinet



Carry Case
CU-SP1-A-BAG-3010



Wall Bracket



Smart Pad Storage



Smart CPR Detection