

CLASSIC-MED BARI

The Classic-Med cushion is designed for comfort and pressure reduction for heavier patients. The castellated foam cells conform to the body redistributing pressure across the surface area. Ideally recommended for those At Risk to High Risk of pressure ulcer development, including light weight patients.



A soft multi stretch, water-resistant cover provides protection from fluid ingress. Welded seams and a protective zip flap support infection control requirements. Easy to wipe down or remove and launder.

FEATURES		BENEFITS	
Bariatric Cushion for Patients up to Hight Risk		At Risk to High Risk	
Comfort for Light Weight Patients		Soft Density Foam Redistributes Pressure from Bony Prominence	
Castellated Foam		Allows Further Immersion and Contour to Body Shape	
Cost Effective Solution		Affordable Solution for Pressure Management	
Multi-Stretch Soft Cover		Soft on the Skin Conforms with the Foam	
High Frequency Welded Seams and Water-Resistant Cover		Reduce Risk of Fluid Ingress, Easy to Wipe Clean	
Rear Zip with Welded Protective Flap		Removable Cover, Launder up to 95 Degrees	
Acti-Fresh Biocide Anti-Bacterial Agent		Integral Anti-Bacterial Protection	
FOAM SPECIFICATION			
Foam Component		Design	Single Castellated Foam
Nominal Density	kg/m3	Nominal Density	N
PRODUCT INFORMATION			
Code	APH072	RISK LEVEL	HIGH RISK
Cushion Size (Sizes available on request)	Made to Order	Weight Capacity	318kg
Warranty	1 year	Product Weight	2.2kg
Safety Standards	BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5		