



Curafix® H

Wide Adhesive Dressing Retention Tape

Product benefits

- Elastic
- Permeable to air and water vapour

Indications

All-over retention of wound dressings in difficult-to dress areas of the body, such as the joints or neck
Retention of drains, catheters etc.

Product Composition

Nonwoven: white, 100% polyester, skin-friendly polyacrylate adhesive (free of colophony and colophony derivatives)



Curafix® H

individually packed, in a folding box

Size (cm x m roll)	L&R REF	Boxes / Cartons
5 x 10	138941	1 / 48
10 x 10	138942	1 / 24
15 x 10	138943	1 / 24
20 x 10	138944	1 / 16
30 x 10	138945	1 / 16