

Curapor® Surgical Non-Woven Island Dressing

Product benefits

- Breathable
- Skin-friendly
- Easy to apply
- Reduced risk of infection (uninterrupted adhesive along all edges)
- Rounded edges for longer adhesion
- Protects against contamination and mechanical influences

Indications

- For the treatment of postoperative, acute and superficial wounds and incisions and lacerations
- Lightly to moderately exuding wounds

Product Composition

Nonwoven base material: white, 100% polyester, polyacrylate adhesive (free of colophony and colophony derivatives)

Wound pad: absorbent layer made of viscose and polyethylene/polypropylene, wound contact layer made from polyethylene

Cover: silicone paper, white



Curapor®
sterile, individually sealed

Size (cm)	Wound pad (cm)	L&R REF	Boxes / Cartons
7 x 5	2.5 x 4	32912	100 / 600
10 x 8	4 x 6	32913	50 / 300
10 x 15	5 x 11	32914	50 / 300
10 x 20	5 x 16	32915	50 / 450
10 x 25	5 x 20	32916	50 / 250
10 x 34	5 x 29	32918	50 / 250