


# danumed® GastroTube



ENSwivel

## Gastrostomy probe with retention balloon for long-term enteral nutrition

Improved product features of the danumed® GastroTube for ease of use and patient comfort:

- **Innovative rotating ring ENSwivel®**, for easy connection and release of ENFit® connections. No twisting of tubes, as only the swivel is turned.
- **New 90° danuPlate** - soft external retention plate with 90° angulation. **Very easy handling!**



## COMFORT

- **New! ENSwivel®** on all ENFit® connectors
- **New! 90° danuPlate** - soft adjustable retention plate for secure and discrete fixation
- Structure of the 90° danuPlate allows excellent **air circulation** in the stoma site area
- **Extensive selection of sizes** from 12 Fr to 24 Fr for the individual and perfect fit
- **Highest patient comfort** due to soft, flexible silicone
- **Length marking** in cm along the tube
- **Display of the balloon volume** on the shaft
- **Convenient accessories** with every GastroTube: a pack of **lubricant** for easy insertion and a **danumed® balloon syringe 10ml** to inflate the retention balloon



## SAFETY

- Gastric replacement feeding tube made of **high-quality medical grade silicone**
- All danumed® GastroTubes are **free of Latex, PVC, DEHP and BPA**
- **Flat rounded tip** allows atraumatic insertion
- **Recessed tip** to prevent trauma
- **ENFit®** connectors according to DIN EN ISO 80369-3 avoid tubing misconnections with other systems
- Available with **ENFit®** and **Y-ENFit®** connectors
- **Radiopaque stripe** along the shaft

## danumed® GastroTubes with ENFit® and Y-ENFit®



REF	Product description	Size in Fr	Pieces/sales unit	Latex-free	PVC-free	DEHP-free	Bisphenol-A-free
DA 1044 0112	danumed® GastroTube with ENFit® and ENSwivel®	12	1 piece	●	●	●	●
DA 1044 0114	danumed® GastroTube with ENFit® and ENSwivel®	14	1 piece	●	●	●	●
DA 1044 0116	danumed® GastroTube with ENFit® and ENSwivel®	16	1 piece	●	●	●	●
DA 1044 0118	danumed® GastroTube with ENFit® and ENSwivel®	18	1 piece	●	●	●	●
DA 1044 0120	danumed® GastroTube with ENFit® and ENSwivel®	20	1 piece	●	●	●	●
DA 1044 0122	danumed® GastroTube with ENFit® and ENSwivel®	22	1 piece	●	●	●	●
DA 1044 0124	danumed® GastroTube with ENFit® and ENSwivel®	24	1 piece	●	●	●	●
DA 1044 0212	danumed® GastroTube with Y-ENFit® and ENSwivel®	12	1 piece	●	●	●	●
DA 1044 0214	danumed® GastroTube with Y-ENFit® and ENSwivel®	14	1 piece	●	●	●	●
DA 1044 0216	danumed® GastroTube with Y-ENFit® and ENSwivel®	16	1 piece	●	●	●	●
DA 1044 0218	danumed® GastroTube with Y-ENFit® and ENSwivel®	18	1 piece	●	●	●	●
DA 1044 0220	danumed® GastroTube with Y-ENFit® and ENSwivel®	20	1 piece	●	●	●	●
DA 1044 0222	danumed® GastroTube with Y-ENFit® and ENSwivel®	22	1 piece	●	●	●	●
DA 1044 0224	danumed® GastroTube with Y-ENFit® and ENSwivel®	24	1 piece	●	●	●	●

Each danumed® GastroTube is individually sterile packed. The danumed® GastroTube is delivered in a sterile peel-pouch and additionally with lubricant and a sterile danumed® balloon syringe 10ml in each single box. 10 boxes in a shipping carton.



INNOVATORS IN ENTERAL CARE

For more information, please contact us on (02) 9882 3666 or at [amsl@amsl.com.au](mailto:amsl@amsl.com.au)

[amsl.com.au](http://amsl.com.au)

AMSL is a subsidiary of Dexcom. PR-100-583 June 2022

