

Diffuser Cushion

Relieve pressure
whilst seated



Memory Foam
2.5 cm

High Density Inner Core
of Traditional Foam
3cm

Diffuser Cushion is designed to:

- Lower pressure and shear forces.
- Provide stability & comfort on chairs.
- Have an absorbent & breathable surface.
- Be easily cleaned and sterilised.
- Ideal as a wheelchair cushion or regular seating.

The diffuser cushion comes in 3 sizes:

- Small • Medium • Large
- Custom Size (with extra cost)

As well as two cover options:


Dura-Fab



Steri-Plus (waterproof)



100% AUSTRALIAN
MADE MEMORY &
TRADITIONAL FOAM



EVEN DISTRIBUTION
OF PRESSURE



IDEAL FOR WHEELCHAIRS,
OUTDOOR AND INDOOR,
FOR OFFICE & CAR



RECOMMENDED BY
AUSTRALIAN HEALTH
PROFESSIONALS