

4H400E ELECTRIC SHOWER TROLLEY

The 4H400E Electric Shower Trolley supports efficient and safe showering and hygiene care of immobile patients. Electrically adjustable positioning promotes correct working posture, while side frames and a soft padded liner provide comfort and support for the patient.



Electric Adjustable Height: Supports ergonomic working positions for carers, through a 400mm range controlled by the hand control.

Stainless Steel Frame: Durable construction enhances product lifespan and reliability.

Soft Shower Liner: Provides patient comfort during showering and care.

Heavy Duty Castors: Three total-lock and one directional 125mm castors for controlled movement and stability.

Removable Battery: Allows fast battery change for continuous operation.



1 YEAR WARRANTY
Components, Hand Control & Cables

5 YEAR WARRANTY
Frame

3 YEAR WARRANTY
Electronics



Max Total Load
(Patient & Water)

Code	Description	Total Length	Total Width	Total Height
BTS147100	4H400E Electric Shower Trolley	1890mm	730mm	780 - 1180mm

ABFY23.11 210325