



REHASENSE



# EXPLORER

## OUTDOOR ROLLATOR

The large outdoor version within the Rehasense rollator range. Weighing 9 kg, with a wide frame and large soft wheels.

# EXPLORER

## OUTDOOR ROLLATOR



REHASENSE

### FOLDABLE AND EASY TO CARRY

**Explorer** folds easily and stays that way with the locking system that doubles up as an ergonomically shaped carrying handle.

### STEADY AND DURABLE

Frame and seat are tested and approved for maximum user weight of 200 kg. The brake mechanism is light yet positive. Soft wheels constructed with double PU layer.

### SIZES

The **Explorer** is available in 2 sizes – Medium and Large – with a range of settings within these sizes



Push handles are height adjustable



Available in 2 sizes – medium and large



Soft wheels constructed with double PU layer

### TECHNICAL INFORMATION

#### Colours:

Brown, black

#### Max. user weight:

200 kg

#### Seat height and article no.:

Large 62 cm seat height

ERBOL600 (Brown)

ERBKL600 (Black)

#### Medium 53 cm seat height

ERBOM550 (Brown)

#### Measurements:

(width x depth x height)

#### Large:

69 x 77 x 91 - 112 cm

Push handle heights: 91 - 112 cm

For body height: 170 - 210 cm

#### Medium:

69 x 77 x 78 - 90 cm

Push handle heights 78 - 90 cm

For body height: 1,50 - 1,80 m

#### Seat width:

53 cm

#### Transportation

measurements folded:

(depth x height x width)

Large : 77 x 91 x 25 cm

Medium : 77 x 84 x 25 cm

#### Weight:

Large : 9,2 kg (+ basket 330 g)

Medium : 9,0 kg (+ basket 330 g)

#### Shopping bag:

Width x depth x height:

33 x 17 x 20 cm

#### Max. load:

5 kg

#### Frame material:

aluminium

#### Wheels:

Front: 285 x 44 mm black soft PU tyre

Rear: 235 x 44 mm black soft PU tyre

### HEIGHT ADJUSTABLE PUSH HANDLES

The **Explorers** push handles are height adjustable, from 78 cm to 112 cm.

The seat heights are 53 and 62 cm with an extra wide 53 cm seat base.

User comfort is ensured with soft wheels.

### HAND GRIPS

Ergonomic shaped hand grips adjustable to the position of the hand.

Smooth operating hand brakes.

### ACCESSORIES

The **Explorer** can be fitted with most of the broad range of Rehasense accessories.



Stays firmly closed with the lock and is easy to open with the push button



Practical easy-to-remove shopping bag. With special designed easy-to-go clips



Maintenance free brakes and TPE-covered unbreakable wheels

**Explorer** is CE-marked and fulfils the requirements of EU directive 93/42 EEC for medical devices.

Copyright © Rehasense® 2019

Rehasense sp. z o.o.  
Sulejowska 45g  
97-300 Piotrków Trybunalski  
Poland

Rehasense Benelux  
Markeweg 81  
8398 GM Blesdijke  
Netherlands

Rehasense Deutschland  
Briener Strasse 25  
47533 Kleve  
Deutschland