

# FLOCARE® FEEDING SETS & TWO PACK CONNECTOR

All Flocare® Feeding Sets are ENFit™ compatible and designed to safely and securely connect enteral tube feeds to feeding tubes and accessories.

## FEATURES

- Flocare® pack feeding sets connect to enteral feed via **ENPlus cross-spike end**: for increased patient safety.
- Flocare® mobile feeding sets are shorter (1580mm) than standard feeding sets (2300mm): for patient convenience and to support mobile feeding.
- The Flocare® two pack connector is a short Y-shaped connector with two ENPlus cross-spike ends: enables connection to two packs of enteral tube feed and/or water at the same time for patient convenience. Quick release clamp may be used to regulate flow.
- All Flocare® feeding sets are compatible with Flocare® containers (40441, 40455, 70139).
- Medication port compatible with ENFit™ syringes: for increased patient safety.
- Made from Latex-free and DEHP-free materials.


### STANDARD FEEDING SETS FOR FLOCARE® INFINITY™ PUMPS

For use with ENPlus Enteral Packs & OpTri Bottles	For use with Enteral Bottles
<b>Flocare® Infinity™ Pack Set – Y Port</b>  Product Code: <b>95347</b> Units Per Carton: 30	<b>Flocare® Infinity™ Bottle Set – Y Port &amp; Drip Chamber</b>  Product Code: <b>101587</b> Units Per Carton: 30
<b>Flocare® Infinity™ Pack Set – No Port</b>  Product Code: <b>95349</b> Units Per Carton: 30	<b>Flocare® Infinity™ Bottle Set – No Port &amp; Drip Chamber</b>  Product Code: <b>90086</b> Units Per Carton: 30
<b>Flocare® Infinity™ Pack Set – 3-way Port &amp; Drip Chamber</b>  Product Code: <b>86514</b> Units Per Carton: 30	<b>Flocare® Infinity™ Bottle Set – 3-way Port &amp; Drip Chamber</b>  Product Code: <b>90084</b> Units Per Carton: 30

## MOBILE FEEDING SETS FOR FLOCARE® INFINITY™ PUMPS

For use with ENPlus Enteral Packs & OpTri Bottles


### Flocare® Infinity™ Pack Mobile Set – Y Port



Product Code: **101586**

Units Per Carton: 30

### Flocare® Infinity™ Pack Mobile Set – No Port



Product Code: **86484**

Units Per Carton: 30

## GRAVITY FEEDING SETS

Not for use with Enteral Feeding Pumps

### Flocare® Gravity Pack Set – Y Port & Drip Chamber



Product Code: **86460**

Units Per Carton: 30


### Flocare® Gravity Pack & Bottle Set – Y Port & Drip Chamber



Product Code: **89755**

Units Per Carton: 30

## FLOCARE® TWO PACK CONNECTOR



Product Code: **69915**

Units Per Carton: 30

### Directions for Use

- Flocare® feeding sets are intended for 'single use only' and should be changed every 24 hours.
- Flocare® two pack connectors are intended for 'single use only' and should be changed every 24 hours.

### Storage

- Store in a clean, dry environment.
- Protect from dust, rain and direct sunlight.

### Ordering Information

To order contact Nutricia Customer Care **1800 889 480**.

For more information call the **Nutricia Careline 1800 438 500**

**NUTRICIA**  
LIFE-TRANSFORMING NUTRITION

Nutricia Australia Pty Ltd.  
Lvl 4, Building D, 12-24 Talavera Road,  
Macquarie Park, NSW, 2113.  
[www.nutriciamedical.com.au](http://www.nutriciamedical.com.au)

NC4845-Oct20 / NUT1180C