


FLOCARE® ACCESSORIES

The Flocare® accessories range includes feeding connectors, adaptors, containers, and pump accessories to aid enteral feeding.

All Flocare® Feeding Tubes, Flocare® Feeding Sets and ENFit™ enteral syringes are ENFit™ compatible.

ADAPTORS	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
 <p data-bbox="225 1055 536 1088">ENFit™ Transition Adaptor</p>	<p data-bbox="644 831 1023 1025">Designed to be inserted into a non-ENFit™ feeding tube to enable connection to an ENFit™ feeding set. Recommended for single patient use for the lifespan of the feeding tube.</p>	<p data-bbox="1091 913 1182 947">122073</p>	<p data-bbox="1310 913 1347 947">50</p>
 <p data-bbox="181 1435 576 1469">Flocare® Transition Step Connector</p>	<p data-bbox="671 1211 995 1406">Designed to be attached to an ENFit™ feeding set to enable connection to a non-ENFit™ feeding tube. Recommended for 24-hour, single use, single patient only.</p>	<p data-bbox="1094 1294 1182 1328">89828</p>	<p data-bbox="1310 1294 1347 1328">30</p>
 <p data-bbox="197 1816 563 1850">Flocare® ENFit™ Bolus Adaptor</p>	<p data-bbox="644 1592 1023 1787">Designed to be inserted into ENPlus enteral feeding packs to enable connection to an ENFit™ syringe to allow for bolus feeding. Recommended for 24-hour, single tube feed pack use.</p>	<p data-bbox="1094 1673 1182 1706">89740</p>	<p data-bbox="1310 1673 1347 1706">30</p>

To order, please call:

Nutricia Customer Care 1800 889 480

CONTAINERS	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
 <p data-bbox="161 461 443 490">Flocare® Container 0.5L</p>	<p data-bbox="555 275 1038 432">Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Recommended for 24 hour, single use only.</p>	<p data-bbox="1134 338 1209 367">40441</p>	<p data-bbox="1385 338 1422 367">10</p>
 <p data-bbox="177 786 427 815">Flocare® Container 1L</p>	<p data-bbox="555 600 1038 757">Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Recommended for 24 hour, single use only.</p>	<p data-bbox="1129 663 1214 692">40455</p>	<p data-bbox="1385 663 1422 692">10</p>
 <p data-bbox="118 1111 486 1140">Flocare® Top-fill Container 1.3L</p>	<p data-bbox="555 902 1038 1093">Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Top-fill design for ease of use. Recommended for 24 hour, single use only.</p>	<p data-bbox="1134 981 1209 1010">70139</p>	<p data-bbox="1385 981 1422 1010">10</p>

MOBILE ACCESSORIES	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
 <p data-bbox="145 1568 459 1597">Flocare® Infinity GO Frame</p>	<p data-bbox="576 1373 1023 1529">Table-top, lightweight, plastic frame that easily fixes to the Flocare® Infinity pump and will allow enteral pack / bottle, Op Tri bottles or Flocare® containers to be attached, for ease of delivery.</p>	<p data-bbox="1134 1440 1214 1469">43247</p>	<p data-bbox="1394 1440 1406 1469">1</p>
 <p data-bbox="161 1892 443 1921">Flocare® Infinity GO Bag</p>	<p data-bbox="563 1697 1034 1854">Designed to optimise patient mobility and product protection when tube feeding. Capacity to carry both the Flocare® Infinity pump, enteral feed and GO frame at the same time.</p>	<p data-bbox="1134 1760 1214 1789">44497</p>	<p data-bbox="1394 1760 1406 1789">1</p>

For more information, please call:
Nutricia Clinical Care Line 1800 060 051



Nutricia Australia Pty Ltd.
 Lvl 4, Building D, 12-24 Talavera Road,
 Macquarie Park, NSW, 2113.
www.nutriciamedical.com.au

NC4656-Apr20 / NUT1141