



SHARPS Disposal Solutions

Eco-friendly

Our sharps container bases are made from 100% recycled plastic, with lids crafted from virgin polypropylene (PP) for superior stability and functionality, ensuring optimal safety during use and transportation. Designed to be puncture-resistant, they provide reliable protection in any healthcare setting.

Quality Assurance

The GuardOn range of containers are manufactured to conform with International Standard ISO 23907-1:2019.

Safety Features

- Easy assembly. Simply push down lid onto its base.
- One-way safety wings (round models) or safety flap (square models) allows for easy single-handed disposal of sharps, yet prevents the outflow of contents
- Opaque-coloured lids to obscure view of contents
- Carry handles for safe transportation
- Temporary closure and permanent closure mechanism with hermetic seal



Pocket Size

THE POCKET provides convenient ease of use for disposal of personal sharps.

Patented double opening system can be used one-handed.

- Side One: pen needle remover
- Side Two: larger opening for disposal of entire pen device or other sharps item.

The Pocket's innovative design and ergonomic shape easily fits into a glove compartment, car console or bag.



Code	Volume	Dimensions (H x W x L)	Qty per Carton
POCKET04L	0.4ml	16.5 x 12.4 x 6.3cm	70

Compact Sharps Containers (2L, 3L, 5L)

2L Round Sharps Container

This 2-piece sharps container features an easy push-down lid assembly for quick and secure setup. Its round, corner-free design reduces the risk of accidental bumps, while the nestable structure saves valuable storage space. The base is made from 100% recycled material, combining functionality with sustainability.

Cap opening has safety wings to prevent accidental spillage of contents.

- Safety wings allow for easy single-handed disposal of sharps, yet prevent the outflow of contents
- Push down narrow end of cap for single-handed temporary closure
- Push down wider end of cap for permanent closure.



Code	Volume (L)	Height (cm)	Diameter (cm)	Qty per Carton	Compatible Bracket
PBS2LECO	2	21.5	14	75	2901900

Compact Sharps Containers (2L, 3L, 5L)

3L Square Sharps Container

This 2-piece sharps container is designed for convenience and safety, featuring an easy push-down lid assembly, a nestable design for efficient storage, and a base made from 100% recycled material.

The container includes a temporary closure with a flip-top opening, an internal safety flap for single-handed sharps disposal while preventing outflow, and a permanent closure secured with just four clicks. An integrated handle ensures easy transport, and an optional mounting bracket is available for added versatility.



Code	Volume (L)	Dimensions (cm)	Qty per Carton	Compatible Bracket
CSP3LECO	3	H 23.7 x W 15.1 x D 15.1	50	2904100

5L Square Sharps Container

This 2-piece sharps container combines safety, sustainability, and convenience. It features an easy push-down lid assembly, a nestable design to save storage space, and a base made from 100% recycled material.

A flip-top opening provides temporary closure, while the internal safety flap enables single-handed sharps disposal and prevents outflow. The container securely seals with just four clicks for permanent closure and includes an integrated handle for easy transport, with an optional mounting bracket available for added flexibility.



Code	Volume(L)	Dimensions (cm)	Qty per Carton	Compatible Bracket
CSP5LECO	5	H 30 x W 16.7 x D 16.7	34	2904500

Medium Volume Sharps Containers (7L)

7L Round Sharps Container

This 2-piece sharps container offers a practical and sustainable solution, featuring an easy push-down lid assembly and a nestable design to save storage space. The base is made from 100% recycled material, and its round design eliminates corners to reduce the risk of accidental bumps. An optional wall mount bracket is available for added convenience.

Cap opening has safety wings to prevent accidental spillage of contents.

- Safety wings allow for easy single-handed disposal of sharps, yet prevent the outflow of contents
- Push down narrow end of cap for single-handed temporary closure
- Push down wider end of cap for permanent closure.



Code	Volume (L)	Height (cm)	Diameter (cm)	Qty per Carton	Compatible Bracket
PBS7LECO	7	19.5	Ø 25	40	2900900

7L Square Sharps Container

This 2-piece sharps container is designed for safety, convenience, and sustainability. It features an easy push-down lid assembly, a nestable design to save storage space, and a base made from 100% recycled material.

The flip-top opening allows for temporary closure, while an internal safety flap enables single-handed sharps disposal and prevents outflow.

The container securely seals with just four clicks for permanent closure and includes an integrated handle for easy transport. An optional mounting bracket is available for added versatility.



Code	Volume (L)	Dimensions (cm)	Qty per Carton	Compatible Bracket
CSP7LBECO	7	H 40 x W 16.7 x D 16.7	27	2904500

Large Volume Sharps Containers (12L, 22L)

12L Round Sharps Container

This 2-piece sharps container offers a practical and sustainable solution, featuring an easy push-down lid assembly and a nestable design to save storage space. The base is made from 100% recycled material, and its round design eliminates corners to reduce the risk of accidental bumps. An optional wall mount bracket is available for added convenience.

Cap opening has safety wings to prevent accidental spillage of contents.

- Safety wings allow for easy single-handed disposal of sharps, yet prevent the outflow of contents
- Push down narrow end of cap for single-handed temporary closure
- Push down wider end of cap for permanent closure.



Code	Volume (L)	Height (cm)	Diameter (cm)	Qty per Carton	Compatible Bracket
PBS12LECO	12	31.3	Ø 25	25	2900900

22L Large Sharps Container

This 2-piece sharps container combines safety, functionality, and sustainability. It features an easy push-down lid assembly, a nestable design to save storage space, and a base made from 100% recycled material.

The flip-top opening provides temporary closure. The container secures permanently with just four clicks and includes molded carry handles for easy transport. A universal mounting bracket is also available for added convenience.



Code	Volume (L)	Dimensions (cm)	Qty per Carton	Compatible Bracket
CSP22LECO	22	H 34.5 x W 37.5 x D 27.2	10	2901200

Cytotoxic Sharps Containers

22L Large Cytotoxic Sharps Container *

This 2-piece sharps container is specifically designed for cytotoxic waste, with a purple color for easy identification. It features an easy push-down lid assembly, a nestable design to save storage space, and a base made from 100% recycled material.

The flip-top opening allows for temporary closure. The container secures permanently with just four clicks and includes molded carry handles for easy transport. A universal mounting bracket is also available for added convenience.



Code	Volume (L)	Dimensions (cm)	Qty per Carton	Compatible Bracket
CSP22LECOPL	22	H 34.5 x W 37.5 x D 27.2	10	2901200

* Non stock item. Samples available upon request. Lead time applies.

Brackets & Accessories

Product Code	Description	To Suit Sharps Container	Qty EA / Carton
2901900 	Powder-coated steel bracket, white. - Mount to vertical pole with clamp (included with screws). - Affix to wall with 2 screws (not included).	2L Round: PBS2LECO	12
2904100 	The white, powder-coated steel cage bracket, combined with the adjustable universal bracket, offers exceptional versatility. This design allows the container to be securely mounted on either a vertical or horizontal pole or directly onto a wall, with screws included for convenience.	3L Square: CSP3LECO	16
2900900 	Powder-coated steel bracket with holes, white. Affix to wall with 2 screws (not included)	7L, 12L Round: PBS7LECO, PBS12LECO	10
2904500 	The white, powder-coated steel cage bracket, combined with the adjustable universal bracket, offers exceptional versatility. This design allows the container to be securely mounted on either a vertical or horizontal pole or directly onto a wall, with screws included for convenience.	5L, 7L Square: CSP5LECO CSP7LBECO	10
2901200 	Universal bracket to mount onto a vertical or horizontal pole, or directly onto a wall (screws included)	Large 22L yellow or purple: CSP22LECO, CSP22LECOPL	50

Contact Information



+612 9868 7708



info.au@munglobal.com



Suite 305, 51 Rawson Street,
Epping NSW 2121, Australia



www.munglobal.com.au

