

## HEALTHWEIGH WHEELCHAIR SCALE, DUAL RAMP AND HANDRAIL



**Product Code:** JE0690

**Product Description:**

Scale - Healthweigh Wheelchair Scale with Two Ramps

**Unit:** EA

### Features

- Movement compensation technology
- BMI function
- Tare/Zero/Hold/Print/kg/lb
- Scroll Down/Up
- 5 digit ¾" LCD/Positive FSTN
- Automatic Reweigh
- Updatable firmware
- Ports: RS232 and USB
- Optional Bluetooth
- Folded, pre-assembly
- 230v/120v input, DC 9V output
- Auto power off
- 6 X 1.5AA batteries
- Low battery indicator

### Specifications

Code - JE0690

Product Size - 1250 x 740 x 1250mm

Platform Size - 800 x 780 x 70mm

Max User Weight - 360kg

Accuracy - 0.1kg

Unit Weight - 43.2kg

Warranty - 2 years