



Trust Quality. Reliable Performance. Cost Effective.

Apogee IC intermittent catheters feature smooth eyelets that are designed to aid comfortable insertion and removal. Non lubricated and not made with natural rubber latex.



Smooth eyelets aid comfortable insertion and removal

Apogee IC
Intermittent Catheter

Apogee IC Intermittent Catheters

Comfortable. Trusted quality.

Apogee IC Features

- Not made with natural rubber latex
- Straight and coude tip options available
- Non-lubricated

Flexible Funnel

- The flexible funnel helps promote easier handling
- Guide notch indicates coude tip position



Packaging

- Easy to open package
- Transparent packaging for ease of viewing

Curved packaging – more discrete and more convenient



Prior to use of Apogee Intermittent Catheters, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

For more information call Customer Care on:

1800 880 851

www.hollister.com.au



Straight Tip Catheter



Tiemann Tip Catheter

Male Catheters

Apogee: Intermittent Catheters

Size	Length	Stock Code/ Replaces	Box QTY
6Fr Firm	38cm	10616/1007	30
8Fr Firm	40cm	10816/1018	30
10Fr Firm	40cm	11016/1034	30
12Fr Firm	40cm	11216/1048	30
14Fr Soft	40cm	1061	30
14Fr Firm	40cm	11416/1064	30
14Fr Curved Packaging Firm	40cm	1065	50
16Fr Firm	40cm	11616/1082	30
18Fr Firm	40cm	11816/1096	30
10Fr Tiemann Tip Firm	40cm	11026	30
12Fr Tiemann Tip Firm	40cm	11226	30
14Fr Tiemann Tip Firm	40cm	11426	30
16Fr Tiemann Tip Firm	40cm	11626	30

Female Catheters

Apogee: Intermittent Catheters

Size	Length	Stock Code/ Replaces	Box QTY
8Fr Firm	15cm	10806/1016	30
10Fr Firm	15cm	11006/1033	30
12Fr Firm	15cm	11206/1047	30
14Fr Firm	15cm	11406/1062	30