

ID SOFT Tilt-In-Space Wheelchair

The ID Soft Tilt-In-Space wheelchair offers various adjustable components, providing maximum adaptability for added comfort. The sturdy steel frame provides extra support and stability in any position for added security.

- Multi-adjustable wheelchair
- Can be tilted and reclined
- Infinite position adjustment of the backrest and seat



The ID Soft Tilt-in-Space Wheelchair provides a smoother ride due to its larger wheels, and yet remains easy to store and transport, also providing optimum climbing ability.



ID SOFT

Tilt-In-Space Wheelchair

The Tilt-in-Space ID Soft Wheelchair can be tilted and reclined into various positions to assist with pressure distribution, correct posture and maximum comfort.

Product Features

- Multi-adjustable wheelchair designed to adapt to most needs of the users requirements
- The wheelchair can be tilted and reclined to help the user maintain correct posture, as well as assisting with even pressure distribution
- Offers maximum comfort for those spending longer periods of time in a wheelchair, and who may also need assistance with their posture
- Height adjustable push handles to provide comfort for the attendant
- Infinite position adjustment of the backrest and seat angle which can be easily operated from the push handle and can be quickly adjusted
- Frame design ensures the wheelchair is steady, whatever the backrest and seat angles are
- Padded adjustable headrest with multi-positional ear cushions and elevating leg rests for additional support and pressure relief
- Anti-tip wheels with integral stepper tube to aid use on kerbs
- Adjustable height backrest and two choices of cushion height
- Long, curved adjustable armrests are fully padded for extra comfort and are easily detachable through a push button
- Back rest reclines from 90° (vertical) to 120°
- Three seat widths are available - 39-44cm (15-17") 44-49cm (17-19") and 49-54cm (19-21")
- Multi-adjustable leg rests can be altered in width, height and depth of calf cushion
- Suitable for use as a seat in a motor vehicle ("Crash Tested" to ISO 7176/19)
- Weight capacity 135kg (21st)
- Accessories available
- This wheelchair has been crash tested providing added safety when being transported in a vehicle

CODE	LENGTH (WITH/ WITHOUT LEGRESTS)	HEIGHT	WIDTH	BACK HEIGHT	SEAT WIDTH	SEAT HEIGHT	SEAT DEPTH	WEIGHT	WEIGHT CAPACITY
Mid-Wheel									
IDSOFTMW16AU	82/108cm (32/42.5")	108cm (43")	58cm (23")	49-64cm (19-25")	39-44cm (15-17")	43cm (17")	40-46cm (16-18")	35kg (77lb)	135kg (21st)
IDSOFTMW18AU	82/108cm (32/42.5")	108cm (43")	63cm (25")	49-64cm (19-25")	44-48cm (17-19")	43cm (17")	40-46cm (16-18")	36kg (79lb)	135kg (21st)
IDSOFTMW20AU	82/108cm (32/42.5")	108cm (43")	68cm (27")	49-64cm (19-25")	48-53cm (19-21")	43cm (17")	40-46cm (16-18")	37kg (81.5lb)	135kg (21st)
Self Propel									
IDSOFTSP16AU	82/108cm (32/42.5")	108cm (43")	58cm (23")	49-64cm (29"-25")	39-44cm (15-17")	43cm (17")	40-46cm (16-18")	33.6kg (74lb)	135kg (21st)
IDSOFTSP18AU	82/108cm (32/42.5")	108cm (43")	63cm (25")	49-64cm (29-25")	44-48cm (17-19")	43cm (17")	40-46cm (16-18")	34.5kg (76lb)	135kg (21st)
IDSOFTSP20AU	82/108cm (32/42.5")	108cm (43")	68cm (27")	49-64cm (29-25")	48-53cm (19-21")	43cm (17")	40-46cm (16-18")	35.3kg (78lb)	135kg (21st)

Accessories

Abduction Wedge

- Can be easily fitted to separate legs and maintain correct posture
- Offers optimum user support whilst remaining extremely comfortable
- Helps stabilise the legs in a position
- Easy wipe clean surface
- Padded with height and length adjustment to ensure user's individualized comfort

CODE:

M73601700AU IDSOFT Abduction Wedge



Standard Tray

- Manufactured from clear plastic enabling the wheelchair user to easily see the front corners of their chair when manoeuvring
- All round 5mm lipped edging prevents items from falling off and retains any accidental spills
- Removable tray to allow easy user transfers with patterned top
- Ideal for reading, working, and for use at mealtimes
- Large cut-out area and depth adjustment gives the user plenty of room and allows positioning close to the tray
- Easy wipe clean surface
- Tray measures a generous 58cm (L) x 37.5cm (D) (23" x 15")

*Not compatible with hemiplegic arms

CODE: M70540800AU IDSOFT Standard Tray



Trunk Support

- Height and angle adjustable, ensuring user's individualized comfort
- Padded cushion support provides added comfort for the user
- Secures the torso at desired angle
- Stays open for easy transfers
- Easy wipe clean surface
- Removable covers for easy cleaning

CODES:

M73601602AU ID Trunk Support (Left Side)

M73601601AU ID Trunk Support (Right Side)



Leaflet code: DDVH_005_IDSOFTH_Wchair