



COVIDIEN

positive results for life™

Kangaroo Joey™ Enteral Feeding Pump

FEATURES

- Water and shock-resistant
- Attitude independent
- "Smart" pump set technology
- Intuitive "ATM-style" ease of operation
- Feeding history up to 72 hours
- DEHP-free pump sets
- Anti-free-flow pump sets
- Longer battery life
- Quick-release pole clamp

Change the world of ambulatory feeding and improve patient outcomes with the pump that feeds and hydrates

An innovation from Covidien that provides feeding and programmable flushing in a completely portable device. The first hydrating ambulatory enteral device available, Covidien's Kangaroo Joey™ pump assures clinicians of accurate and reliable function while providing patients with freedom and mobility.

Unique ambulatory feeding and flushing capabilities

- Deliver consistent and reliable enteral nutrition and hydration while patients are on the go.

Fully portable design with backpack

- Allows for convenient, discrete use.

Compact, lightweight construction

- Ideal for use by all ages requiring home or ambulatory care.

Easy, intuitive interface

- Ensures accurate and reliable delivery of preset volumes.

New ambulatory-specific feeding set design

- Assists in ensuring proper set-up and formula administration.



Easy-to-use pump set



Proximal spike Set



1000 mL feed and flush pump set



1000 ml feed only pump set



Kangaroo™ safety screw spike

www.kangarooepump.com



FEATURES

- Pumping Mechanism – Rotary peristaltic
- Operating Range – 1 - 400 mL/hour in 1 mL/hr increments
- Volume to be Delivered – 1 - 3000 mL in 1 mL/hr increments
- Accuracy – +/- 7%
- Occlusion Pressure – 15psi
- Battery Life – 18 hours at 125 mL/hr
- Memory – 72 hours
- Dimensions – 4.1 x 5.1 x 3.6
- Weight – 1.7 lbs., w/pole clamp 2.3 lbs.
- Alarms – 8 options
- Multilingual – 15 Languages
- Tubing length: 9" (from pump to feeding container).

ORDERING INFORMATION

| Code Number | Description | Case Quantity |
|-------------|--|---------------|
| 383400 | Kangaroo Joey Enteral Feeding Pump w/ Pole Clamp | 1 |
| 383493 | Kangaroo Joey Pole Clamp | 1 |
| 383494 | Kangaroo Joey Flexible Pole Clamp | 1 |
| 770024 | Small Kangaroo Joey Backpack | 1 |
| 383491 | Kangaroo Joey Power Cord | 1 |
| FO10506 | Kangaroo Joey Battery | 1 |

PUMP SETS

| | | |
|--------|---|----|
| 764655 | Kangaroo Joey Proximal Spike Set | 30 |
| 762055 | Kangaroo Joey 500 mL Pump Set | 30 |
| 763656 | Kangaroo Joey 1000 mL Pump Set | 30 |
| 764669 | Kangaroo Joey Spike Set w/ 1000 mL Flush Bag | 30 |
| 763662 | Kangaroo Joey 1000 mL Feed Set with 1000 mL Flush Set | 30 |

RECERTIFICATION SET

| | | |
|--------|------------------------------------|---|
| 776150 | 100 mL Burette Recertification Set | 5 |
|--------|------------------------------------|---|

KANGAROO FEEDING SETS FOR USE WITH NESTLÉ®* SPIKERIGHT™* SYSTEM

| | | |
|--------|--|----|
| 765555 | Kangaroo Joey Safety Screw Spike | 30 |
| 765662 | Kangaroo Joey Safety Screw Spike w/1000 mL Flush Bag | 30 |
| 775762 | Safety Screw Spike Adapter Cap | 30 |

©*Nestlé is a registered trademark of the Société des Produits Nestlé S.A., Vevey, Switzerland.
 ™*SpikeRight is a trademark of the Société des Produits Nestlé S.A., Vevey, Switzerland.
 COVIDIEN, COVIDIEN with logo and ™ marked brands are trademarks of Covidien AG or an affiliate.
 © 2009 Covidien. All rights reserved.

H6268-5M-0609