

# THE NEXT STEP IN COMFORT

Our new all-in-one  
JOBST® forMen



DEFY GRAVITY

Our new all-in-one JOBST forMen uniquely combines our best features for men with new, enhanced benefits. It is specially engineered to meet men's daily needs and to provide all-day comfort and relief.

# OUR ALL-IN-ONE PRODUCT FOR MEN, UNIQUELY COMBINES OUR BEST FEATURES



**SoftFit**  
Ensures socks stay up  
for elevated comfort



**Durable fabric**  
Designed for long-lasting  
compression



**StayFresh**  
Effective odour control



**Specially constructed sole**  
For cushioning




**New anatomic toe area**  
For enhanced comfort  
and fit



# HELP YOUR PATIENTS TAKE THE NEXT STEP IN COMFORT.

With a wide range of options, our new JOBST forMen makes it easier to provide men with the all-day comfort and relief they need.

### STYLE



Knee High      Thigh High

### TOE OPTIONS



Closed toe

### COMPRESSION CLASSES\*

15-20 mmHG

20-30 mmHG

30-40 mmHG

\*The mean compression for an average ankle size.

### COLOUR OPTIONS



Black      Khaki


### TOP BANDS



SoftFit (Knee high)      Micro dotted (Thigh high)

## SIZE


Knee High



	S	M	L	XL	LFC	XLFC
cC	28-38 cm	30-42 cm	32-46 cm	34-50 cm	46-61 cm	46-61 cm
cB	18-21 cm	21-25 cm	25-29 cm	29-33 cm	25-30 cm	30-36 cm

Length D: 38 - 45 cm

Thigh High



	S	M	L	XL
cG	40-62 cm	46-70 cm	54-78 cm	60-81 cm
cC	28-38 cm	30-42 cm	32-46 cm	34-50 cm
cB	18-21 cm	21-25 cm	25-29 cm	29-33 cm

Length G: 68.5 - 80 cm

# WE JUST TOOK THE NEXT STEP IN COMFORT. HELP YOUR PATIENTS DO THE SAME.

We created our new JOBST forMen to help more men with venous conditions experience the all-day lightness and comfort they deserve. By combining our best features for men with new enhanced benefits, we engineered a new all-in-one compression sock that gives men the all-day comfort and support they need.

Our new JOBST forMen is soft yet supportive, light yet durable, and distinct yet discreet - one simple recommendation that meets men's daily needs. Recommend new JOBST forMen to your patients.



Length	S	M	L	XL	LFC	XLFC	Colour
	Closed Toe						
Knee High	60523-40	60523-41	60523-42	60523-43	60523-44	60523-45	● Black
Knee High	60524-40	60524-41	60524-42	60524-43	60524-44	60524-45	● Khaki
Thigh High	60423-20	60423-21	60423-22	60423-23			● Black
Thigh High	60423-60	60423-61	60423-62	60423-63			● Khaki
Knee High	60529-40	60529-41	60529-42	60529-43	60529-44	60529-45	● Black
Knee High	60530-40	60530-41	60530-42	60530-43	60530-44	60530-45	● Khaki
Thigh High	60425-20	60425-21	60425-22	60425-23			● Black
Thigh High	60425-60	60425-61	60425-62	60425-63			● Khaki
Knee High	60535-40	60535-41	60535-42	60535-43	60535-44	60535-45	● Black
Knee High	60536-40	60536-41	60536-42	60536-43	60536-44	60536-45	● Khaki
Thigh High	60427-20	60427-21	60427-22	60427-23			● Black
Thigh High	60427-60	60427-61	60427-62	60427-63			● Khaki



JOBST®,  
an Essity brand

Always read the label and follow the directions for use.

ESS1134  
08/2025

#### Essity Health & Medical

PO Box 337, Mount Waverley,  
Victoria, Australia 3149  
medical.essity.com.au  
jobst.com.au

JOBST® Customer Service AU  
T 1300 998 810 F 1300 998 820

JOBST® Customer Service NZ  
T 0508 998 810 F 0508 998 820