

# JOBST® ULCERCARE



JOBST® UlcerCARE™ enables effective management of hard to heal wounds on the lower leg while offering ease of use for patients and care givers - the perfect combination for providing patients with an improved quality of life.

## JOBST® UlcerCARE™ Compression System:

- Provides graduated compression of 40mmHg, the recommended pressure for healing these types of wounds
- This 2-in-1 stocking system is designed for patient self-management to facilitate easy changing of the wound dressing and to allow for improved personal hygiene (e.g. bathing and skin care)
- JOBST® UlcerCARE Liner provides mild compression (15-20mmHg) for the non-ambulatory patient to help manage lower leg oedema during bed rest
- Outer stocking comes with a zipper on the left or right side for comfortable closure on the opposite side of the wound. It is removed at night and applied first thing in the morning

## Sizing & Compression Classes

Size	Ankle		Calf		JOBST® UlcerCARE Zipper Stocking and Liners (1 stocking/2 liners)		
	Zipper on Left Side	Zipper on Right Side	JOBST® UlcerCARE Replacement Liners (3 liners)				
Small	18 - 21 cm	29 - 35 cm	60428-00	60428-20	60428-60		
Medium	21 - 24 cm	34 - 41 cm	60428-01	60428-21	60428-61		
Large	24 - 27 cm	40 - 47 cm	60428-02	60428-22	60428-62		
X-Large	27 - 30 cm	46 - 54 cm	60428-03	60428-23	60428-63		
XX-Large	30 - 33 cm	48 - 56 cm	60428-04	60428-24	60428-64		
XXX-Large	33 - 36 cm	50 - 58 cm	60428-05	60428-25	60428-65		
XXXX-Large	36 - 39 cm	52 - 60 cm	60428-06	60428-26	60428-66		

JOBST UlcerCARE is supplied with one outer stocking and two liners without a silicone top band. Replacement liners are supplied in packs of three.

