

JOBST® soSoft

Circular knit


Ready-To-Wear (RTW)

- JOBST® soSoft compression socks provide graduated compression to help support circulation
- With gentle compression of 8-15mmHg, JOBST® soSoft compression socks can help tired and achy legs
- Suitable for longer periods of sitting and standing
- Provide relief from minor oedema





Sizing & Compression Classes

8-15 mmHg
Supportive Compression



S	M	L	XL	Colour
73679-11	73679-21	73679-31	73673-05	● Black
73680-12	73680-22	73680-32	73673-10	● Sand

Sizing		Ankle cB	Calf cC
	S	18 - 21 cm	28 - 38 cm
	M	20 - 24 cm	30 - 42 cm
	L	22 - 26 cm	32 - 46 cm
	XL	24 - 28 cm	34 - 50 cm

US Sizing Women		
	S	4.5 - 6.5
	M	7 - 9
	L	9.5 - 11
	XL	11.5 - 13

Shoe size alone is not an indicator for correct sizing. Measurements of calf and ankle circumference are mandatory for the proper performance of this product

