



By your side
and beyond

From the Cardinal Health™ Kangaroo™ Portfolio

The future of enteral feeding

The Kangaroo OMNI™ Enteral Feeding Pump provides patients flexibility in their feeding preferences to make their home feeding journey easier.

Tube feeding is an adjustment that can present patients with emotional and physical challenges. Gastrointestinal issues or food allergies may cause discomfort, while feelings of isolation can bring on anxiety. Thick formula can help address some of these challenges, but many feeding pumps on the market aren't able to successfully deliver thick formula.

Thick formula is defined as enteral fluids of smooth consistency that are categorized as Level 2, 3 or 4 Drinks within the International Dysphagia Diet Standardization Initiative (IDDSI) Framework. **The Cardinal Health™ Kangaroo OMNI™ Enteral Feeding Pump is designed to deliver thick formula** to help provide enteral feeding patients more options to meet their unique tube feeding needs. Our Thick Formula Feeding Sets are innovative, including larger feeding set inner tube diameters to accommodate more viscous formulas (refer to Figure 1).




Figure 1

*Tube drawings not to scale

Features:



Capable of delivering thick formula
(Levels 2, 3 and 4 on IDDSI framework)



When used with Thick Formula Feeding Sets, pump has accuracy of +/- 10% or 1.0mL/hr, whichever is larger.



Stores up to 30 days of feeding history




Thick Formula Feeding Sets are not made with DEHP, latex or BPA

IDDSI Framework:

The first internationally recognized organization that methodically developed a framework to categorize the viscosities of various foods and liquids. It provides two testing methods:

- 1. Flow Test:** Uses objective measurement tool to determine liquid thickness (Levels 0-3)
- 2. Fork Drip Test:** Assesses whether moderately and extremely thick drinks and liquidized foods (Levels 3 and 4) flow through fork prongs



© The International Dysphagia Diet Standardization Initiative 2019

With the Kangaroo OMNI™ Enteral Feeding Pump, give your patients greater flexibility to receive the holistic nutrition they need and prefer.

For more information visit
cardinalhealth.com.au/kangaroo

Important Information: Prior to use, refer to the Instructions for Use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings, and precautions.

Cardinal Health Australia 503 Pty Ltd | cardinalhealthanz@cardinalhealth.com | +61 (0) 2 8311 6500 | cardinalhealth.com.au
© 2024 Cardinal Health. All Rights Reserved. CARDINAL HEALTH and the Cardinal Health LOGO are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. MSCR-702_Mar-2024.
Whilst endeavours have been taken to ensure the content of this material is correct at the time of publishing, Cardinal Health Australia 503 Pty Ltd (ABN 22 607 002 460) and its related entities or its agents, do not warrant the accuracy of any information, documents, images or illustrations included in this document and reserve the right to change the content of the document at any time. Interested parties must rely on their own enquiries.

