

# Luer Lock Syringe 1mL

## Low Dead Space

### PRODUCT DATA SHEET

#### DESCRIPTION

m|devices 1mL Luer Lock Syringe is intended to deliver an accurate metered dosage throughout the administration of a hypodermic injection. The Low Dead Space allows for minimal residual volume due to the displacement plunger tip. The high transparency barrel allows the syringe contents to be clearly seen and the finger grip on the plunger end ensures stability when administering contents.


#### FEATURES AND BENEFITS

- The Low Dead Space allows for minimal residual volume due to the displacement plunger tip
- High transparency barrel
- Luer Lock for a secure connection to corresponding hypodermic or other component
- Finger grip on plunger end ensures stability and comfort during injection

#### PRODUCT DATA

- Syringe Volume: 1mL
- Material: Polypropylene barrel and plunger, latex free rubber gasket and silicone lubricant.
- Sterilisation: EO
- DEHP, Latex and PVC Free
- Shelf Life: 5 Years
- ARTG Number: 280799

#### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV189001	Luer Lock Syringe, Low Dead Space	1mL	100