




Suction Tubes and Accessories
Tracheobronchial Suction Catheters Product Guide

Medtronic



CONTENTS

Tracheobronchial Suction Catheters	4
Specimen Traps	8


TRACHEOBRONCHIAL SUCTION CATHETERS

TOUCH-TROL™

Argyle™ bronchial suction catheter for avoiding suction trauma

- Uniformly distributed flow
- Colour coded
- Suitable for aspiration of mucous from the tracheobronchial tract
- Medical grade PVC
- Tube for bronchial suction with bevelled, atraumatic tip
- Available with the Touch-Trol™ suction control system
- With one central hole and two side eyes to protect the mucosa from invagination.

Tip Detail



1180851147

Product Code	Description	Colour	Diameter Fr	Length ¹ cm	Qty/Case
1180851063		Light green	6	44	100
1180851089	Standard suction catheters with vacuum breaker connector and 2 side eyes	Blue	8	50	100
1180851105		Black	10	50	100
1180851121		White	12	50	100
1180851147		Green	14	50	100
1180851162		Orange	16	50	100
1180851188		Red	18	50	100

AERO-FLO™

Argyle™ bronchial suction catheter with atraumatic tip

- Available with Touch-Trol™ suction control system and funnel connector
- Slide-Tex™ external finish
- Medical grade PVC
- Tube for bronchial suction with atraumatic "ring" tip
- With one central hole and four side holes to protect the mucosa from invagination and avoid suction trauma
- Uniformly distributed flow

Tip Detail



8888250142

Product Code	Description	Colour	Diameter Fr	Length ¹ cm	Qty/Case
8888250092		Lavander	5	25	50
8888250118		Blue	8	40	50
8888250126	Atraumatic suction catheters with Touch-Trol™ Vacuum breaker connector and integral ring tip at distal end	Black	10	52	50
8888250134		White	12	52	50
8888250142		Green	14	52	50
1180251744		White	12	51	50
1180251751	with coude tip	Green	14	51	50

STANDARD SUCTION CATHETER

- Touch-Trol™ suction control system
- Medical grade PVC
- Tube for bronchial suction with atraumatic “ring” tip
- With one central hole and four side holes to protect the mucosa from invagination and avoid suction trauma
- Uniformly distributed flow




LOOP PACK

Product Code	Description	Size	Length	Qty/ Pack
PHS07	Suction catheter standard	7FG	39cm	100
PHS08	Suction catheter standard	8FG	39cm	100
PHS10	Suction catheter standard	10FG	47cm	100
PHS10L	Suction catheter standard – Loop pack	10FG	56cm	50
PHS12	Suction catheter standard	12FG	47cm	100
PHS12L	Suction catheter standard – Loop pack	12FG	56cm	50
PHS14	Suction catheter standard	14FG	47cm	100
PHS14L	Suction catheter standard – Loop pack	14FG	56cm	50
PHS16	Suction catheter standard	16FG	47cm	100
31208	Suction catheter standard	8FG	39cm	100
31110	Suction catheter standard	10FG	56cm	100
31112	Suction catheter standard	12FG	56cm	100
31114	Suction catheter standard	14FG	56cm	100

SUCTION CATHETER WITH DEPTH MARKINGS

- Depth marking (cm)
- Touch-Trol™ suction control system
- Medical grade PVC
- Tube for bronchial suction with atraumatic "ring" tip
- With one central hole and four side holes to protect the mucosa from invagination and avoid suction trauma
- Uniformly distributed flow



DEPTH MARKINGS (CM)




Product Code	Description	Size	Length	Qty/Pack
31305	Suction catheter standard with depth markings	5FG	39cm	100
31306	Suction catheter standard with depth markings	6FG	39cm	100
31307	Suction catheter standard with depth markings	7FG	39cm	100
31308	Suction catheter standard with depth markings	8FG	39cm	100
31310	Suction catheter standard with depth markings	10FG	56cm	100

SPECIMEN TRAPS

LUKENS

Argyle™ System for suctioning bronchial secretions and collecting mucous

- Whistle-tip suction catheter with Slide-Tex™ finish
- Transparent 20 cc specimen trap graduated at 2 cc intervals
- Removable container cap
- Additional cap
- Suction tube with universal connector
- Sticky label for patient identification
- Preassembled



8888258608

Product Code	Description	Diameter Fr	Length ¹ cm	Qty/Case
8888258608	Mucus trap container with DeLee catheter	14	20	50

SPECIMEN TRAP

Argyle™ specimen trap for bronchial secretions

- Also useable with a bronchoscope
- Preassembled
- Screw-on cap
- Removable male connector
- Transparent specimen trap graduated 1 cc intervals
- Sticky, pressure-sensitive label for patient identification



8884724500

Product Code	Description	Capacity cc	Qty/Case
8884724504	To be used with bronchoscope	70	50
8884724500	DeLee Neo safe bubble	40	50

DeLEE


Argyle™ System for suctioning bronchial secretions and collecting mucous

- Suction catheter with Slide-Tex™ finish
- Transparent specimen trap graduated at 2 cc intervals
- Removable container cap
- Tube with moulded mouthpiece
- Preassembled
- Medical grade PVC



8888257527

Product Code	Description	Diameter Fr	Capacity cc	Qty/Case
8888257527	Mucus trap container with DeLee catheter	10	20	50



¹ as referenced on the packaging which is effective length or straight part of catheter.

Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road
Macquarie Park NSW 2113
Australia
Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 668 670

medtronic.com.au
medtronic.co.nz

Medtronic New Zealand Ltd
Level 3 - Building 5, Central Park Corporate Centre
666 Great South Road,
Penrose, Auckland 1051
New Zealand
Fax: +64 9 918 3742
Toll Free: 0800 377 807

© Medtronic.
2017 All Rights Reserved.
1192-01-18 # 4632-012018