

# Mobile Shower Commode



Corrosion resistant accessories

Universal design for easy fitment

4x lockable castors

Height adjustable

Height adjustable leg supports

All aluminium frame with powder coated finish

Swing back arm supports

Padded PU seating and comfort components

Pan and lid

Swing away leg supports

Flip up foot rests

This Mobile Shower Commode is designed for aged care, community and home care environments. The frame and its major components are constructed from corrosion resistant aluminium and powder coated for extra protection. All accessories, including pins, nuts and bolts are constructed with corrosion resistant materials to help ensure the strength and longevity of the product when used as intended.

## Features

- Multi-function design for improved usability
- Universal fitment suits most toilet types
- Height adjustable for improved ergonomics and fitment
- Padded accessories for user comfort
- Corrosion resistant construction materials
- 4x lockable castors for improved safety and stability
- Swing back arm supports for easy transfer
- Enclosed pan and lid with improved spill protection
- Flip up foot rests with removable leg supports
- Flat pack design for compact and easy storage and/or transport

## Specifications

PRODUCT CODE	BL0030
COLOUR	White
PRODUCT WEIGHT	12 kg
OVERALL HEIGHT	930-1005 mm
OVERALL WIDTH	550 mm
OVERALL DEPTH	840 mm
SEAT HEIGHT	495-570 mm
SEAT WIDTH	440 mm
SEAT DEPTH	435 mm
WIDTH BETWEEN ARMRESTS	450 mm
CASTORS	4 x Lockable
FRAME	Aluminium
SAFE WORKING LOAD	136 kg
CONFORMS TO	AS/NZS 3973:2009
ARTG	176831
WARRANTY	1 Year