



Morpheus Air Static Hybrid™ MAS Replacement Mattress

Product Code - APH235

Specifications

Dimensions: 198cm x 88cm x 15.6cm

Maximum Weight Limit: 39 stone / 247.7 Kg

Risk Factor: Very High Risk

Warranty: 5 years cover / 5 years foam

Mattress Weight: 16.5 Kg

Safety Standards: BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5



Apollo Healthcare Technologies Limited, Holme Street, Liversedge, West Yorkshire WF15 6JF
Tel: +44 (0) 1924 614567 Fax: +44 (0) 1924 607480 Email: sales@apollo-ht.co.uk www.apollo-ht.co.uk

www.apollo-ht.co.uk

Morpheus Air Static Hybrid™ MAS Replacement Mattress

Product Code / Dimensions:
APH235 / 198cm x 88cm x 15.6cm



The **Morpheus Air Static Hybrid™** mattress is designed for users considered to be at **“Very High Risk”** of pressure ulcer development. It combines the benefits of effective air displacement technology with modern visco-elastic foam for excellent pressure reduction.

With high levels of comfort offered and requiring no power unit/pump to operate the system, it is ideal for community and palliative care settings, whilst the clinical benefits and ease of use make it suitable in an acute environment. Incorporating 10 air and foam filled cells, a unique valve system allows Reactive Air to flow through the system to displace and adjust to the patient's body weight and movement. It is equipped with a visco-elastic foam top layer.

- All component parts are interchangeable and replaceable reducing environmental impact and offering a very cost effective solution
- Vapour permeable, multi way stretch fabric cover high frequency welded seams and deep flap zip protection
- Designed for users at very high risk, or have existing tissue damage
- Fire retardant cover to Standards BS7175 Crib 5
- The mattress cover can be wipe down or completely removed and laundered up to 95°C to aid decontamination protocols

