


The Morpheus Air Zero is an Advance 4 Zoned Mattress Support System offering both Full Dynamic and Semi Dynamic Therapy.


- Offering therapeutic immersion that simulates dynamic therapy whilst unpowered.
- The Addition of a Pump converts the mattress to dynamic mode when higher therapy is required



- The Morpheus Air Zero is suitable for all patient risks, offering reactive immersion therapy by means of innovative self-adjusting cell technology.
- The Morpheus Air Zero can be converted to dynamic therapy by the addition of a pump when a higher mode of therapy is required.
- The Morpheus Air Zero behaves like a Dynamic System when unpowered, making it the mattress of choice for smokers when powered devices cannot be used.


Innovative Foam Topper works in conjunction with the Top Cover designed to address shear and friction increasing comfort and support.

Multi Stretch, Vapour permeable cover.




Pressure Mapping

Powered



Non-Powered



Product Code	APH166Z
Risk Category	Very High
Max User Weight	250 kg
Cell Cycle	1 in 2 (10 Minutes)
Innovative Foam Topper	✓
Foam in TPU Cell	✓
Width	900mm
Depth	150mm
Mattress Weight	14 kg
Weight of Pump	1.7 kg
Fault Alarm	✓
Power Failure Alarm	✓
Easy Grip Handles	✓
Cable Management	✓
CPR	✓
Fully Profiling Flexible Side Walls	✓
4 Zoned Heel, Sacral, Torso and Head	✓
Available as a Bariatric Option	✓
Multi Stretch Vapour Permeable Cover	✓

Cell ZONE 1
HEEL Section

Cell ZONE 2
SACRAL Section

Cell ZONE 3
TORSO Section

Cell ZONE 4
HEAD Section

