

GloveON™

Cyton Long Cuff Exam Gloves

Specifically developed for use in cytotoxic and chemotherapy environments, GloveON Cyton gloves, offer reliable barrier protection without compromising comfort or dexterity. Designed with healthcare professionals in mind, they combine strength, durability, and ease of use to provide the highest confidence in high-risk settings.

Features:

- Purple Colour – easily identifiable for cytotoxic use, reducing cross-contamination risks
- High Strength Nitrile – excellent puncture and chemical resistance.
- Extended Cuff – added forearm protection in high-risk procedures.

Suitable for:

- Handling cytotoxic and chemotherapy drugs
- Oncology and chemotherapy wards
- Pharmacy compounding
- Hospital and day oncology units



Features		
<p>Distinctive Colour</p>	<p>Protection & Compliance</p>	<p>Comfort & Dexterity</p>

Product Code	Size	Sterile	Colour	Pcs/ Unit	Unit/ Ctn
CYTON LONG CUFF GLOVES					
CYL811XS	XSmall	Non-Sterile	Indigo	100	10
CYL811SS	Small	Non-Sterile	Indigo	100	10
CYL811MM	Medium	Non-Sterile	Indigo	100	10
CYL811LL	Large	Non-Sterile	Indigo	100	10
CYL811XL	XLarge	Non-Sterile	Indigo	100	10
CYL811XXL	XXLarge	Non-Sterile	Indigo	80	10