

# Purus 23L Touch Series II

## - B class Autoclave

Purus recognizes the importance of quality, reliability, and value for money.

Our range of autoclaves brings easy-to-use functionality, contemporary styling, technical quality, and the reassurance of warranty protection and service support.

This means that your business can ensure the highest quality of sterilization standards which you and your clients deserve.

### 23 Litre B class

- LCD display window displays sterilisation process
- 4 button control with smart touch keys
- High precision, temperature sensors precise to 0.1 degree, pressure sensor accurate to .001bar
- Multiple record output, print USB or SD card
- Stand by setting
- Custom cycle function
- Hand lock door with pressure lock system
- Smart maintenance system
- Water quality test function
- European standard, with 3 times pre and drying vacuum
- Latest design steam generator
- Sterilizing Temp 121 °C 134 °C
- Sterilizes:
  - Packed/unpacked solid or hollow items, porous items
  - Removes air from hollow area of any instrument
- Australian standard AS/NZ 4815-2006, CE
- **TGA ARTG Registration No 222444 Class 11B**
- Total compliance for your sterilization needs



23 Litre B class Touch Autoclave

### Product Code

purus\_AC-23Btouch Series 11

### Product Description

#### Chamber Size

- Diameter - 245mm
- Depth - 470mm

#### External dimensions

- 67cm D x 47cm W x 40cm H

#### Weight

- 55kg.

# 23L Touch Series II - B class Autoclave

## Installation

- Simply plug in and use – self-installed by customer
- Bowie Dick, Vacuum and Helix test Installed
- Accessories supplied:
  - USB connection
  - chamber rack
  - 3 trays
  - door spanner
  - tray tool
  - spare door seal
  - level

## Service

- Autoclave Calibration / Validation to be completed annually by a Purus Authorised Service Technician
- Quality and cost effective service supported by warranty
- Repair and Maintenance agreements available
- No hidden upkeep service requirements


## Warranty

- 3 year back to base warranty on Autoclave

## Finance

We have a panel of finance providers (including Medfin, Westpac Health, FinLease and Investec) and will introduce you to them for a competitive financing package.

### Autoclave calibration graph



### Contact Details

Contact 1300 Purus 9  
PO Box 2111 SUNBURY VIC 3429  
p: +61 3 8358 5557  
e: info@purushealth.com.au  
www.purushealth.com.au

## Call Purus Now

- To discuss your Autoclave needs
- For best pricing