



REHASENSE



SERVER ^W

WITH GREAT STYLE

The **Server W** is based on the stylish Router rollator, including all the same features and functions. The soft and comfortable arm supports and the brake handles are multi adjustable in depth, angle and height.

SERVER^W

WITH GREAT STYLE



REHASENSE

FOLDABLE

It is foldable, staying firmly closed making it is easy to transport.

STEADY AND DURABLE

The frame and seat are tested and approved on users with a maximum weight of 150 kg. Light and reliable brakes act on the TPE covered unbreakable wheels. The optimal and rigid design of the frame and the fore arm support ensures the best possible usage comfort (lightness, turning circle and steadiness).

SIZES

The **Server W** is available in size large.

HEIGHT ADJUSTABLE

The **Server W** is easily height adjustable in the range from 90 – 115 cm.

ACCESSORIES

Server W can be fitted with most of the Rehasense accessories.




Ergonomic shaped hand grips adjustable to the position of the hand. Smooth operating hand brakes.



Seat height 62 cm.




Maintenance free brakes and TPE-covered unbreakable wheels.



Stays firmly closed with the lock and is easy to open with the push button.



Practical easy-to-remove shopping bag with special designed easy-to-go clips.



Unique triangle profile ensures the durability up to 150 kg.

TECHNICAL INFORMATION

Standard colours:



Wine red Sober grey

Max. user weight:

150 kg

Measurements:

(width x depth x height)

55 x 66 x 90 – 115 cm

Seat height:

62 cm

Weight:

10,8 kg (+ basket 300 g)

Frame material:

Aluminium

Shopping bag:

(width x depth x height)

32 x 15 x 28 cm

Item no.:

NSWRL60UK (wine red)

NRSGL60UK (sober grey)

Server W is CE-marked and fulfils the requirements of EU directive 93/42 EEC for medical devices.

Copyright © Rehasense® 2019

Rehasense UK
Mobility House
Marsh Lane Industrial Park
Hayle, Cornwall
TR27 5JR

www.rehasense.com

info@rehasense.com