



COMFORT & SUPPORT

PATIENT LIFTING SLINGS





COMFORT AND SUPPORT

Solari Medical specialises in providing quality products to the aged and health care industries of Australia.

Manufactured to Australian Standards our extensive product range aims to provide the highest level of comfort and support.



GENERAL PURPOSE SLING

- Divided legs can be crossed in the centre, uncrossed or positioned under thighs.
- Suitable for patients with good head control.
- Available in solid polyester or mesh.

SIZE	A	B	C	D	POLY	MESH
SMALL	800MM	1150MM	550MM	400MM	SGPP1	SGPM1
MEDIUM	1000MM	1320MM	600MM	480MM	SGPP2	SGPM2
LARGE	1060MM	1410MM	610MM	560MM	SGPP3	SGPM3



GENERAL PURPOSE SLING WITH HEAD SUPPORT


- For patients needing full body and neck support.
- Divided legs can be crossed in the centre, uncrossed or positioned under thighs.
- Available in solid polyester or mesh.

SIZE	A	B	C	POLY	MESH
SMALL	1200MM	550MM	700MM	SGHP1	SGHM1
MEDIUM	1500MM	650MM	1030MM	SGHP2	SGHM2
LARGE	1650MM	700MM	1130MM	SGHP3	SGHM3

FOR CUSTOM SLINGS AND SERVICING CONTACT SOLARI MEDICAL




SURGICAL HOUSE
 Medical & Surgical Supply Specialists
surgicalhouse.com.au
surgical@surgicalhouse.com.au
 (08) 9381 4199



DELUXE COMMODE SLING

- Soft padding throughout sling.
- Excellent back support and comfort.
- Suitable for patients who generally have good trunk and hip control.
- Long divided legs.
- Available in polyester solid only.

SIZE	A	B	C	D	POLY
SMALL	770MM	980MM	670MM	360MM	SCDP1
MEDIUM	870MM	1010MM	620MM	430MM	SCDP2
LARGE	1000MM	1140MM	750MM	580MM	SCDP3
X-LARGE	1100MM	1190MM	800MM	630MM	SCDP4



TOILET/STANDING SLING

- Soft padding throughout sling.
- Only use with a standing lifter.
- Allows for safe and comfortable chair to toilet transfer.
- Suitable for patients who can weight bear and who have appropriate upper body strength.
- Allows patient to be easily brought to a standing position.
- Has a security belt with buckle for extra safety.
- Available in polyester solid only.

SIZE	A	B	C	POLY
SMALL/MED	910MM	380MM	360MM	STSP1
LARGE/XL	1050MM	400MM	375MM	STSP2

**FOR CUSTOM SLINGS AND SERVICING
CONTACT SOLARI MEDICAL**



SURGICAL HOUSE
Medical & Surgical Supply Specialists
surgicalhouse.com.au
surgical@surgicalhouse.com.au
(08) 9381 4199


COMFORT AND SUPPORT



HAMMOCK SLING

- Non-slip backing prevents sling from sliding under patient.
- Can be left under patient, reducing manual handling.
- Excellent back and neck support.
- Available in solid polyester only.

SIZE	A	B	C	POLY
MEDIUM	480MM	1200MM	860MM	HSCP2
LARGE	580MM	1250MM	960MM	HSCP3
X-LARGE	680MM	1300MM	1060MM	HSCP4



FULL BODY SLING

- Excellent back support and comfort.
- For patients who need full body and neck support.
- Available in solid polyester or mesh.

SIZE	A	B	C	POLY	MESH
SMALL	1250MM	400MM	800MM	SFBP1	SRFM1
MEDIUM	1370MM	450MM	830MM	SRFP2	SRFM2
LARGE	1450MM	550MM	900MM	SRFP3	SRFM3

**FOR CUSTOM SLINGS AND SERVICING
CONTACT SOLARI MEDICAL**



SURGICAL HOUSE
Medical & Surgical Supply Specialists
surgicalhouse.com.au
surgical@surgicalhouse.com.au
(08) 9381 4199

COMFORT AND SUPPORT



FEATURES

STRONG BONDING WAIST STRAP

- Anti-Snag & Anti-Linting Touch Tape
- Stronger Bonding
- Less Abrasive




KEYHOLE ATTACHMENT

HDPE keyhole attachment available.



COLOUR-CODED LIFTING STRAPS

High strength 50mm seatbelt webbing.



SERVICE & MAINTENANCE LABEL

To ensure you can have an accurate service history. As per Australian Standard.

2019												
2020												
2021												
2022												
	Jan	Feb	Mar	Apr	May	Jun	Jul	Aug	Sep	Oct	Nov	Dec

6 MONTH SERVICING SCHEDULE

MADE TO AUSTRALIAN STANDARDS

AS 10535-2011

SAFE WORKING LOAD 300KG



	P 08 9336 7469 F 08 9443 3034 E enquiries@solari.com.au W www.solari.com.au
COMFORT AND SUPPORT	
SWL 300KG MANUFACTURED TO AUSTRALIAN STANDARDS AS 10535-2011 ARTG NO. 226474	
CARE INSTRUCTIONS MACHINE WASH TO 60° - COLD RINSE - DRIP DRY OR TUMBLE DRY TO MAX 40° INSPECT SLING REGULARLY FOR WEAR AND TEAR PRIOR TO USE DO NOT USE IF ANY FAULTS ARE FOUND PLEASE REFER TO USER INSTRUCTIONS IF IN DOUBT	
SLING TYPE	
<input type="checkbox"/> TOILETING/STANDING	<input type="checkbox"/> COMMODE DELUXE/HYGIENE ACCESS
<input type="checkbox"/> GENERAL PURPOSE	<input type="checkbox"/> GP WITH HEAD SUPPORT
<input type="checkbox"/> COMFORT	<input type="checkbox"/> FULL BODY
<input type="checkbox"/> HAMMOCK WITH NON-SLIP BACKING	<input type="checkbox"/> CUSTOM MADE
SLING ATTACHMENT	
<input type="checkbox"/> KEY HOLE	TO BE USED WITH KEYHOLE YOKE HOIST
<input type="checkbox"/> WEBBING LOOP	TO BE USED WITH STANDARD YOKE HOIST
SLING SIZE	
<input type="checkbox"/> S RED	<input type="checkbox"/> M BLUE
<input type="checkbox"/> L GREEN	<input type="checkbox"/> XL BLACK
<input type="checkbox"/> POLY NAVY	<input type="checkbox"/> MESH NAVY
DATE OF MANUFACTURE	
MONTH	1 2 3 4 5 6 7 8 9 10 11 12
YEAR	2018 2019 2020
BATCH #	18 A B C D
COMMENCEMENT DATE:	

WOVEN LABELS

The Care Instructions and Service & Maintenance labels are woven, fade resistant and made to Australian Standards.

**FOR CUSTOM SLINGS AND SERVICING
CONTACT SOLARI MEDICAL**



SURGICAL HOUSE
Medical & Surgical Supply Specialists
surgicalhouse.com.au
surgical@surgicalhouse.com.au
(08) 9381 4199

COMFORT AND SUPPORT