

Solvaline® N

Low Adherent Dressing

Product benefits

Highly absorbent viscose/polyester core
Polyethylene film perforated on both sides for rapid absorption of exudate

Low-adherent: allows virtually painless and atraumatic dressing changes

Film coating on both the top and underside allows use on either side

Can be cut to size

Indications

Exuding wounds

Wounds in the post-operative phase

First aid, e.g., of cuts and abrasions

2nd degree burns

Product composition

Absorbent core: 80% viscose, 20% polyester

Film: polyethylene



Solvaline® N
sterile, individually sealed

Size (cm)	L&R REF	Boxes / Cartons
5 x 5	280004	100 / 2700
10 x 10	280005	100 / 1200
10 x 20	280006	50 / 600
20 x 30	280007	25 / 300
5 x 5	280008	25 / 2700
10 x 10	280009	10 / 720