



SURGICAL HOUSE
Healthcare Consumables & Equipment

STANDARD CONTOUR CUSHION

The Standard Contour Cushion (4" cell height) is Star Cushion Products premier cushion for the prevention and treatment of pressure ulcers. Flexible, interconnected 4" air cells allow for deep immersion into the cushion providing constant low pressure across the contact seating surface. In addition, the cell moves and shapes to the contact area providing a low shear, low friction surface ideal for the prevention and healing of up to and including a Stage IV pressure ulcer.



HCPCS CODE - E2624

Dual Valve



HCPCS CODE - E2624

Four Valve



4" cell height with cover

- Air adjustable 4" cells offer deep immersion to provide a high level of pressure reduction for tissue protection and treatment.
- Available in dual and four valve/compartments to effectively address seating challenges such as pelvic obliquities and sacral sitting, to enhance posture control and to improve overall function.
- Light Weight, Approximately 3 pounds./1.36 kilograms.
- Available in Pediatric and Adult Sizes.
- 24 Month Manufacturers Limited Warranty.
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene rubber.
- No Weight Limit.



4" cell height
HCPCS CODE - E2622