


Suction Catheters Mully with Vacuum Control



Area of application

Suction catheters are used for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via a tube before being connected to vacuum source.

The products are used for single use only.

Product range

Description	CH	Connector color	Effective Length mm	Pack.	Standard (SH A 76) Reference	SOFT (SH A 64) Reference
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 490mm						
Mully Metric with vacuum control, straight tip tip, CH6	6	Light green	490	1/100/600	70000218	
Mully Metric with vacuum control, straight tip tip, CH8	8	Light blue	490	1/100/600	70000418	
Mully Metric with vacuum control, straight tip tip, CH10	10	Black	490	1/100/600	70000518	70400118
Mully Metric with vacuum control, straight tip tip, CH12	12	White	490	1/100/600	70000618	70400218
Mully Metric with vacuum control, straight tip tip, CH14	14	Green	490	1/100/600	70000718	70400318
Mully Metric with vacuum control, straight tip tip, CH16	16	Orange	490	1/100/600	70000818	70400418
Mully Metric with vacuum control, straight tip tip, CH18	18	Red	490	1/100/600		70400518
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 570mm						
Mully Metric with vacuum control, straight tip tip, CH10	10	Black	570	1/100/600	70000928	70400628
Mully Metric with vacuum control, straight tip tip, CH12	12	White	570	1/100/600	70001028	70400728
Mully Metric with vacuum control, straight tip tip, CH14	14	Green	570	1/100/600	70001128	70400828
Mully Metric with vacuum control, straight tip tip, CH16	16	Orange	570	1/100/600	70001228	70400928
Mully Metric with vacuum control, straight tip tip, CH18	18	Red	570	1/100/600	70001328	70401028
Suction Catheters Mully Metric with vacuum connector, Curved tip, Effective length 570mm						
Mully Metric with vacuum control, straight tip tip, CH6	6	Light green	570	1/100/600	90000128	
Mully Metric with vacuum control, straight tip tip, CH8	8	Light blue	570	1/100/600	90000228	
Mully Metric with vacuum control, straight tip tip, CH10	10	Black	570	1/100/600	90000328	
Mully Metric with vacuum control, straight tip tip, CH12	12	White	570	1/100/600	90000428	
Mully Metric with vacuum control, straight tip tip, CH14	14	Green	570	1/100/600	90000528	
Mully Metric with vacuum control, straight tip tip, CH16	16	Orange	570	1/100/600	90000628	
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 330mm						
Mully Metric with vacuum control, straight tip tip, CH8	8	Light blue	330	1/100/600	12500428	
Mully Metric with vacuum control, straight tip tip, CH10	10	Black	330	1/100/600	12500528	
Mully Metric with vacuum control, straight tip tip, CH12	12	White	330	1/100/600	12500628	
Mully Metric with vacuum control, straight tip tip, CH14	14	Green	330	1/100/600	12500728	
Mully Metric with vacuum control, straight tip tip, CH16	16	Orange	330	1/100/600	12500828	

Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 5 years

Suction Catheters Mully with Vacuum Control



Certifications	
CE-marking	CE 0123
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485, Medical Device Single Audit Program (MDSAP)
European Medical Device Nomenclature (EMDN)	A06010103 CANNULAS AND TIPS FOR APIRATION AND IRRIGATION, SINGLE-USE
Global Medical Device Nomenclature (GMDN)	10749

Instructions for use	
Included?	Yes

Material specifications	
Contains latex?	No
Contains DEHP?	No
Contains cotton?	No
Contains paper?	Packaging only
Tube	Polyvinyl chloride
Connector	Styrene-butadiene copolymers

Packaging information									
REF	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions, mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m ³		
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 490mm									
70000218	1	630x40	100	640x175x120	600	660x375x370	0.092	8.0	9,000
70000418	1	630x40	100	640x175x120	600	660x375x370	0.092	8.5	9,000
70000518	1	630x40	100	640x175x120	600	660x375x370	0.092	9.0	9,000
70000618	1	630x40	100	640x175x120	600	660x375x370	0.092	9.5	9,000
70000718	1	630x40	100	640x175x120	600	660x375x370	0.092	10.0	9,000
70000818	1	630x40	100	640x175x120	600	660x375x370	0.092	11.5	9,000
70400118	1	630x40	100	640x175x120	600	660x375x370	0.092	9.0	9,000
70400218	1	630x40	100	640x175x120	600	660x375x370	0.092	9.5	9,000
70400318	1	630x40	100	640x175x120	600	660x375x370	0.092	10.0	9,000
70400418	1	630x40	100	640x175x120	600	660x375x370	0.092	11.5	9,000
70400518	1	630x40	100	640x175x120	600	660x375x370	0.092	12.0	9,000
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 570mm									
70000928	1	713x40	100	720x175x120	600	740x375x370	0.103	9.5	9,000
70001028	1	713x40	100	720x175x120	600	740x375x370	0.103	10.5	9,000
70001128	1	713x40	100	720x175x120	600	740x375x370	0.103	11.3	9,000
70001228	1	713x40	100	720x175x120	600	740x375x370	0.103	12.8	9,000
70001328	1	713x40	100	720x175x120	600	740x375x370	0.103	13.7	9,000
70400628	1	713x40	100	720x175x120	600	740x375x370	0.103	9.5	9,000
70400728	1	713x40	100	720x175x120	600	740x375x370	0.103	10.5	9,000
70400828	1	713x40	100	720x175x120	600	740x375x370	0.103	11.3	9,000
70400928	1	713x40	100	720x175x120	600	740x375x370	0.103	12.8	9,000
70401028	1	713x40	100	720x175x120	600	740x375x370	0.103	13.7	9,000

Suction Catheters Mully with Vacuum Control



Packaging information continued									
REF	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions, mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m ³		
Suction Catheters Mully Metric with vacuum connector, Curved tip, Effective length 570mm									
90000128	1	713x40	100	720x175x120	600	740x375x370	0.103	9.5	9,000
90000228	1	713x40	100	720x175x120	600	740x375x370	0.103	10.5	9,000
90000328	1	713x40	100	720x175x120	600	740x375x370	0.103	11.3	9,000
90000428	1	713x40	100	720x175x120	600	740x375x370	0.103	12.8	9,000
90000528	1	713x40	100	720x175x120	600	740x375x370	0.103	9.5	9,000
90000628	1	713x40	100	720x175x120	600	740x375x370	0.103	10.5	9,000
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 330mm									
12500428	1	485x33.5	100	500x175x120	600	520x375x370	0.072	7.9	12,000
12500528	1	485x33.5	100	500x175x120	600	520x375x370	0.072	8.3	12,000
12500628	1	485x33.5	100	500x175x120	600	520x375x370	0.072	8.5	12,000
12500728	1	485x33.5	100	500x175x120	600	520x375x370	0.072	9.3	12,000
12500828	1	485x33.5	100	500x175x120	600	520x375x370	0.072	9.7	12,000

GTIN Codes			
REF	Primary package	Secondary package	Tertiary package
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 490mm			
70000218	06936881154089	16936881154086	26936881154083
70000418	06936881154096	16936881154093	26936881154090
70000518	06936881154102	16936881154109	26936881154106
70000618	06936881154119	16936881154116	26936881154113
70000718	06936881154126	16936881154123	26936881154120
70000818	06936881154133	16936881154130	26936881154137
70400118	06936881154140	16936881154147	26936881154144
70400218	06936881154157	16936881154154	26936881154151
70400318	06936881154164	16936881154161	26936881154168
70400418	06936881154171	16936881154178	26936881154175
70400518	06936881154188	16936881154185	26936881154182
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 570mm			
70000928	6936881159596	16936881159593	26936881159590
70001028	6936881159619	16936881159616	26936881159613
70001128	6936881159626	16936881159623	26936881159620
70001228	6936881159633	16936881159630	26936881159637
70001328	6936881159640	16936881159647	26936881159644
70400628	6936881159657	16936881159654	26936881159651
70400728	6936881159664	16936881159661	26936881159668
70400828	6936881159671	16936881159678	26936881159675
70400928	6936881159688	16936881159685	26936881159682
70401028	6936881159695	16936881159692	26936881159699

Suction Catheters Mully with Vacuum Control



GTIN Codes cont			
REF	Primary package	Secondary package	Tertiary package
Suction Catheters Mully Metric with vacuum connector, Curved tip, Effective length 570mm			
90000128	6936881159794	16936881159791	26936881159798
90000228	6936881159800	16936881159807	26936881159804
90000328	6936881159817	16936881159814	26936881159811
90000428	6936881159824	16936881159821	26936881159828
90000528	6936881159831	16936881159838	26936881159835
90000628	6936881159848	16936881159845	26936881159842
Suction Catheters Mully Metric with vacuum connector, Straight tip, Effective length 330mm			
12500428	6936881160097	16936881160094	26936881160091
12500528	6936881160110	16936881160117	26936881160114
12500628	6936881160127	16936881160124	26936881160121
12500728	6936881160134	16936881160131	26936881160138
12500828	6936881160141	16936881160148	26936881160145

Packaging material	
Primary package	Coated paper, polypropylene/polyethylene film
Secondary package	Single-wall corrugated fiberboard
Tertiary package	Double-wall corrugated fiberboard

Suction Catheters Mully with Vacuum Control



Storage and handling	
	THIS WAY UP
	Keep away from rain
	Fragile; handle with care
	Keep away from sunlight
	Do not resterilize
	Do not re-use
	Do not use if package is damaged and consult electronic instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or use electronic instructions for use

SPONSORED IN AUSTRALIA BY

Uniqmedical Pty Ltd
PO Box 1117
Windsor VIC 3181

1300 55 80 79
info@uniqmedical.au

DISTRIBUTED BY

UNI Q MEDICAL GMBH
Schützengraben 20
8200 Schaffhausen
Switzerland
info@uniqmedical.net
www.uniqmedical.net



Well Lead Medical Co., Ltd.
C-4 Jinhu Industrial Estate
Hualong, Panyu Guangzhou
511434 P.R. China



Shanghai Intl. Holding
Corp. GmbH
Eiffestrasse 80
20537 Hamburg
Germany



QS Engineering AG
Erlenstrasse 31
CH - 4106 Therwil
Switzerland

**UK
Responsible
Person
(UKRP)**

SUNGO Certification
Company Limited
3rd Floor, 70 Gracechurch
Street, London, EC3V 274HR
England, United Kingdom