


tg® grip Tubular Support Bandage



Product benefits

- supple material with uniform compression effect
- air-permeable
- light padding effect
- application with virtually no creases
- knitted seamlessly
- widthwise stretch
- skin-friendly due to high percentage of cotton
- washable up to 60°C and sterilisable by steam (134°C)
- compression can be enhanced with a second layer of material or by using a smaller size

Indications

- for support and retention with a light compression effect

Product composition

- 85% cotton, elastodiene (latex); covered with a polyamide yarn
- contains latex



To note

- washable up to 60°C
- sterilisable by steam (134°C)
- compression can be enhanced with a second layer of material or by using a smaller size

tg® grip 10 m roll in a cardboard dispenser

Size	REF	Shipping Units (pcs.)
A (4.6) for children's hands, children's forearms	24 120	1/24
B (6) for children's arms, children's lower legs	24 121	1/16
C (6.5) for (slim) arms, children's legs	24 122	1/16
D (7.5) for medium-sized arms, lower legs	24 123	1/16
E (8.75) for large arms, legs, slim thighs	24 124	1/12
F (10) for legs, medium-sized thighs	24 125	1/12
G (12) for legs, for large thighs	24 126	1/8
J (17.5) for (small) trunks	24 127	1/6
K (21.5) for (medium-sized) trunks	24 128	1/4
L (32.5) for (large) trunks	24 129	1/4