



Suction Tubing (Sterile)

Our sterile Suction Tubing is available in flexible, ribbed and special soft versions and comes in a variety of lengths. Each tube comes with standard hexangular connectors attached.

This product is well suited for surgical use.

- ⊗ Flexible, ribbed and special soft tubing options
- ⊗ Choice of tubing lengths
- ⊗ Hexangular connects
- ⊗ Kink and crush resistant tubing that maintain lumen patency
- ⊗ Securely packaged with double wrap
- ⊗ Single use, sterile and ready to use

Independent testing exceeds minimum requirements for suction tubing degree of collapse as per AS 2120.3 Clauses 5.2.2



Suction Tubing			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
AN050000	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	2m	80
AN050001	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	3m	50
AN050002	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	4.5m	40
AN050003	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	6m	35

Table continued on next page

Suction Tubing (Sterile)

Suction Tubing			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
AN050004	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	8m	25
AN050009	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	1m	50
AN051000	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	2m	80
AN051001	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	3m	50
AN051002	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	4.5m	40
AN051003	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	6m	35
AN051004	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	8m	25
AN054002	Suction Tubing - Ribbed, ID6.8mm OD10mm, 2 connectors	2m	80
AN054006	Suction Tubing - Ribbed, ID6.8mm OD10mm, 2 connectors	4m	50
AN056001	Suction Tubing - Special Soft, ID6mm OD10mm, 1 connector	3m	55
AN056006	Suction Tubing - Special Soft, ID6mm OD10mm, 1 connector	4m	50
AN057001	Suction Tubing - Special Soft, ID6mm OD10mm, 2 connectors	3m	55
AN057006	Suction Tubing - Special Soft, ID6mm OD10mm, 2 connectors	4m	50
AN059000	Suction Tubing - Special Soft, Stepped Connector	0.48m	120