

AirSpiral™ Tube and Auto-fill Chamber Kit

REF MYAIRVOKIT1

for myAirvo™ Series humidifiers



Introducing the MYAIRVOKIT1, AirSpiral Tube and Auto-fill Chamber Kit.

The MYAIRVOKIT1, comprising of an AirSpiral tube and a water chamber, is for delivery of humidified respiratory gases. These are to be used with the myAirvo Series humidifiers in homes and long-term care facilities.

AirSpiral Tube

The AirSpiral tube utilizes an advanced two-spiral design to protect against condensate.



Auto-fill Chamber

To be used with a water bag e.g. 900PT401



Specifications

Product code	MYAIRVOKIT1
Product description	myAirvo AirSpiral Tube and Auto-fill Chamber Kit
Quantity	Each
Usage	Maximum 60 days; Single patient multiple use
Humidified High Flow System	myAirvo Series humidifiers
Compatible Patient Interfaces	Optiflow™ + Series Interfaces 2-Pack (MYOPT9SMALL, MYOPT9MEDIUM, MYOPT9LARGE, MYOPT9TRACHE, MYOPT9MASK) Optiflow Junior 2 Home Nasal Cannulas (OJR416HM, OJR418HM) Optiflow + Series Interfaces 1-Pack (OPT942E, OPT944E, OPT946E, OPT970E, OPT980E) Optiflow Junior Nasal Cannulas (OPT316, OPT318)

Available in the kit below:



MYAIRVOKIT1

Find out more at

www.fphcare.com/myairvo

Some products may not be available in your country. Please contact your Fisher and Paykel Healthcare representative.